

# Lucitone® Intensive Colors

## Quick Start Guide



Grind the denture to roughen the surface at the location where the characterization will be done.


**Note:** Blow with oil-free, compressed air to clean the surface.



Add small portion of the desired Lucitone® Intensive Color(s) powder to the cold cure acrylic powder and stir, then add the cold cure denture base liquid and mix thoroughly. The material will become thicker as it sets on the bench.



Wet the prepared surface with the cold curing liquid.



When appropriate consistency is reached, apply the intensive colored material to the prepared denture.



Cure the characterized denture in a pressure pot according to the DFU of the cold cure denture base material.

- For Lucitone® HIPA Denture Base, cure by placing the characterized denture in a pressure pot with water at  $113^{\circ}\text{F} \pm 2^{\circ}\text{F}$  ( $45^{\circ}\text{C} \pm 1^{\circ}\text{C}$ ) and immediately apply a minimum of 20 psi (1.4 bars) pressure. Cure for 30 minutes.



Repeat process as needed to achieve desired characterization.



Finish and polish as normal.

# Lucitone® Intensive Colors

## Experience Excellence with High-Pigment Colors

Lucitone® Intensive Colors provide esthetic excellence with versatile characterization options. Intensive Colors are highly pigmented acrylic powders that enable the technician to mix with self-curing acrylic powder or liquid systems, and individually characterize a denture. The Lucitone® Intensive Color Kit features 5 base colors that can be used to create a variety of shade options.

- + Mixing Tray
- + Spatula
- + Five Color Shades




- + Yellow
- + Brown
- + Red
- + White
- + Blue



## Ordering Information

Lucitone® Intensive Colors		Size	Shade
905990	Lucitone® Intensive Color Kit	Kit: 5 colors, mixing tray, spatula	5 color kit
905992	Lucitone® IC Powder	10 grams	Red
905993	Lucitone® IC Powder	10 grams	Yellow
905994	Lucitone® IC Powder	10 grams	White
905995	Lucitone® IC Powder	10 grams	Brown
905996	Lucitone® IC Powder	10 grams	Blue

 Manufacturer:  
Dentsply Sirona Prosthetics  
570 West College Avenue, York, PA 17401  
1-800-243-1942  
prosthetics.dentsply.com

DeguDent GmbH  
Rodenbacher Chaussee 4  
63457 Hanau-Wolfgang  
Germany  
+49 6181 59-50