

CLINICAL & CDT CATALOGUE

SALES FREEPHONE
0800 009 2777

Supplying the dental profession since 1960

whw
Behind every smile

TERMS OF TRADING

GENERAL

- 1. In these terms of trading “the Company” means WHW Plastics Ltd.
- 2. Our previous terms of trading are hereby cancelled.
- 3. The placing of an order shall be conclusively deemed an acceptance of these terms and any purported variation or exclusion shall be of no effect unless accepted by the company in writing.

PRICES

- 4. All prices shall be subject to alteration without notice and goods will be sold at the Company's price current on date of despatch.
- 5. All prices shall be deemed to be exclusive of Value Added Tax at the appropriate rate from time to time in force, postage and/or carriage and/or delivery charges (if applicable), insurance and all applicable taxes, duties, tariffs and charges of any nature of whatsoever imposed in any country or territory either directly or indirectly in respect of the sale or supply of goods or payment for them.

PAYMENT

- 6. Payment of all invoices, unless marked otherwise, shall be made by the last day of the month following month of invoice, e.g. purchases during January due for payment 28th February.
- 7. In the event of any default in making payment as aforesaid and without prejudice to the Company's other rights, the company shall be entitled to charge interest on any outstanding balance at the rate of 2% per month above the Base Rate of Natwest Bank from time to time in force and for the period commencing with the day upon payment should have been made as aforesaid and ending on the day which payment is actually made.

EARLY SETTLEMENT DISCOUNT

- 8. Any early settlement discount will be clearly marked on invoices and must be paid within 14 days of invoice.

GOODS RETURNED FOR CREDIT

- 9. Prior consent of the Company must be obtained before any goods will be considered for credit. Subject to this consent, full credit (excluding carriage costs) will be given providing all of the following conditions are met:-
 - (a) The invoice/Delivery Note number and date are quoted when application is made for the Company's consent.
 - (b) The goods are returned carriage paid by the Customer or, with the prior consent of the Company , via a duly authorised representative of the Company.
 - (c) The goods are returned within14 days of delivery to the Customer.
 - (d) The goods are in the same condition as they were in at the time of despatch by the. Company to the Customer. Failure to comply with any of the above conditions will result in a partial credit or no credit at all at the sole discretion of the Company.

SPECIAL ORDER ITEMS

- 10. Special Order items are defined as those not appearing in our catalogue or any amendment thereto or products that are manufactured to the customers specifications. Special Order items may not be returned for credit under any circumstances (unless faulty) nor may Special Orders be cancelled once placed.

NON-DELIVERY

- 11. The Company shall have no liability in respect of the goods lost or damaged in transit unless the Customer notes such shortage loss or damage on the delivery note and gives written notice thereof to the Company and the Carrier with 7 days or such shorter period as may be required by the Carriers Conditions of Carriage. The Company's liability in respect of the goods lost or damaged in transit shall be limited to repairing or replacing such goods.

OWNERSHIP/TITLE

- 12. No property in or title to any goods shall pass to the Customer until the Customer has satisfied all monies owing to the Company by the Customer on whatever account.

RISK/INSURANCE

- 13. Risk passes to the Customer on delivery/collection of the goods. The Customer shall be responsible for and shall indemnify the Company against loss or damage of any sort until such time as he becomes the owner of goods. This applies equally where goods have been supplied on loan or approval.

QUALITY OR FITNESS

- 14. All items, conditions and warranties (whether implied or made expressly) whether by the Company or its servants or agents or otherwise (other than those warranties set out in the current edition of the Company's catalogue) relating to the quality and/or fitness for the purpose of the goods or any of the goods excluded.

MANUFACTURERS INSTRUCTIONS

- 15. The Company shall have no liability in respect of any losses arising of the Customers failure to strictly and fully comply with any manufacturers instructions which may accompany any items detailed in our Catalogue or otherwise supplied by the Company.

05 ACCESSORIES

10 ACID ETCHANTS

08 AMALGAM SEPARATION

10 ANAESTHETICS

11 DENTURE BASE, RELINE & REPAIR

16 HYGIENE

17 IMPRESSION

24 ORTHODONTIC

26 PREVENTIVE

30 PROVISIONAL


32 RESTORATIVE

36 CENTRIFORM

41 CENTRI AFTER CARE

42 PACKAGING

43 CENTRI



AWARD-WINNING PRODUCTS

GELATO PROPHYLAXIS PASTE

2013, 2014 & 2015 Dental Advisor Best Propyl Cup
Page 28

GORILLA SEAL™

Dental Advisor 4 + Rated
Page 31

FASTRAY™ LC

Dental Advisor 4.5 + Rated
Page 31

AEGIS R PIT & FISSURE SEALANT WITH ACP

Dental Advisor 4 + Rated
Page 29

AEGIS R CROWN & BRIDGE CEMENT WITH ACP

Dental Advisor 4 + Rated
Page 32

ACCESSORIES

☐ **Coolley Pic™**
Double End Pic

☐ **Evacuators**
High Volume Oral Evacuators

☐ **E-Z Squeeze™**
Universal Dispenser

☐ **Flexo®**
Saliva Ejectors

☐ **Masque®**
Silicone Lubricant and Separator

☐ **“My Tooth” Boxes®**
Storage Boxes

☐ **Novelty Keychains**
Keychains

☐ **Pressure Pot**
Hydraulic Water Press

☐ **Syringes**
5 ml Disposable Syringes

☐ **TC Retractors**
Tongue and Cheek Retractors

☐ **Visual Systems**
Bookkeeping Forms

☐ **Wonder Wedges®**
Interproximal Wooden Wedges

GENERAL DENTAL NECESSITIES

#0921329	Applicator Brushes 2.5 in. (50)
#0921116	Applicator Brushes (100)
#0921117	Handle (1)
#0921336	Wax-Coated Mixing Cups (80 ml, 3 oz) (100)
#0921785	Plastic Mixing Cups (30 ml) (100)
#0921540	All-Purpose Mixing Pad (6.375 in. x 4.25 in.) (39 sheets)
#0921563	Composite Mixing Pads (2.75 in x 2.125 in) (2 pads)

#0921564	Double-Ended Plastic Spatulas (100)
#0921089	3 ml Pipette Eyedrop-Type Dispenser (15)
#0921867	5 ml Disposable Plastic Syringes (50)
#0921545	Static Mixing Tips 1:1 (2 in.) (100)
#0921546	Static Mixing Tips 4:1 (1 in.) (100)
#0921547	Static Mixing Tips 1:1 (1.25 in.) (100)
#0921548	Static Mixing Tips 1:1 (1 in.) (100)

COOLEY PIC™

Coolley Pic is an accessory for removing temporary restorations, copper bands, inlays, bridges, root tips, and any operation where a strong, pointed instrument is required. Made of stainless steel.

#0921317 Double End Pic



EVACUATORS

Disposable, high volume non-vented oral evacuators protect soft tissues for maximum patient comfort while suctioning blood, saliva and debris. Fits all standard vacuum lines. Latex-free.

Recommended for use with: TC Retractors

#16651 White (100)

#16652 Blue (100)

#16653 Green (100)





E-Z SQUEEZE™ UNIVERSAL DISPENSER & STATIC MIXING TIPS

An economical instrument for dispensing all cartridge-bound impression materials.

Recommended for use with: DentuSil, Ultra Trim, SuperSil and Re'cord



#0921809 E-Z Squeeze Universal Dispenser
#0921313 1:1 EZ Squeeze Dispenser
#0921312 1:1/2:1 EZ Squeeze Dispenser

#0921809 E-Z Universal Dispenser fits the following:



MIXING TIPS

#0921414 1:1 Static Mixing Tips for 4.5 in cartridges (48)
#0921699 4:1 E-Z Squeeze Static Mixing Tips (15)





FLEXPLO®

Flexos are flexible hygienic saliva ejectors that will adapt to any mouth. They maintain their shape and do not collapse when bent. The specially designed soft, smooth tips provide patient comfort and will not pull tissue. Their non-clog design and consistent shape retention will ensure a perfect free flowing fit. Available in clear and four bright, translucent colours: blue, green, red, and yellow. Clear Flexo are available with either a fixed or removable tip.

#1662213	Rainbow Pack (1000) 200 ea. fixed tips	
#1662214	Clear with Clear removable tips (100)	case: 10
#1662215	Clear with White fixed tips (100)	case: 10

#1662216	Blue with White fixed tips(100)	case: 10
#1662217	Yellow with White fixed tips (100)	case: 10
#1662218	Green with White fixed tips (100)	case: 10
#1662219	Red with White fixed tips (100)	case: 10



MASQUE®

Masque is a silicone lubricant and coating agent that is insoluble in the mouth and keeps acrylic from bonding to the teeth. Masque is unaffected by monomers and heat. It cannot run, is tasteless, odourless, non-oily, and completely compatible with oral tissue. Masque can be applied whenever a lubricant or coating agent is required. Use wherever you do not want acrylic to bond. To remove Masque, simply wipe clean and rinse.

#0921530	1 oz (29 g) Syringe with Plunger
----------	----------------------------------



“MY TOOTH” BOXES®

“My Tooth” Boxes are clear-lidded pink or blue boxes specially designed to hold teeth. These economical, brightly coloured boxes make a great gift from the tooth fairy.

Dimensions: 1.125 in x 1.125 in x 0.875 (when closed)

#1662220P	Pink Boxes (100)
#1662220B	Blue Boxes (100)



NOVELTY KEYCHAINS

Stainless steel novelty keychains available in 3 styles. Perfect as gifts for your staff and patients!

#0921793	Mini Perforated Impression Tray
#0921794	Mini Extraction Forceps
#0921795	Mini Dental Mirror



PRESSURE POT

Bosworth's Pressure Pot is a hydraulic water press used in creating bubble-free relines, repairs, splints, temporaries, and orthodontic appliances. It guarantees a dense, super smooth, porosity-free surface on any self-curing acrylic with just 10-15 p.s.i. Create a temporary with Pressure Pot to increase strength by reducing voids and porosity. Pressure Pot does not require electricity, heat, or air pressure. It is lightweight and easily transportable. Comes with a DVD.

Recommended for use with: Masque™

Total weight of unit: 6 lbs 2 oz • Inside Diameter: 5 in/12.7 cm • Outside Diameter: 5.625 in/14.3 cm
Inside Height: 3.5 in/8.9 cm • Overall Height: 7.5 in/19.1 cm • Depth: 5.5 in/14 cm • Width: 6 in/15.2 cm

Presse à eau hydraulique • Hydraulische Wasserpresse • Pressa idraulica ad acqua • Prensa hidráulica de água • Prensa de agua hidráulica

#092135	Pressure Pot
---------	--------------



WATCH THE
TECHNIQUE VIDEO!

SYRINGES

Disposable plastic syringes (5 ml) specially designed for use with rubber base and silicone impression materials. The barrel is translucent for viewing exact amount of material loaded into the syringe. The syringes are perfectly sized for easy, proper placement of material in the patient's mouth. The bendable tip allows for easier, more accurate placement of materials. Syringes are easy to load, hygienic, and save time on clean up.

Einwegspritze (5 ml) • Jeringuillas Disponibles (5 ml) • Une Seringue à Usage Unique (5 ml)

#0921867	Syringes (50)
----------	---------------



TC RETRACTORS

The perfect complimentary item to the high-volume oral evacuator is the TC Retractor. The disposable soft tip retractor fits securely onto the evacuator to retract the tongue and/or cheek. The built in “filter” prevents inlays, crowns and bridges from being lost into the evacuation system.

Recommended for use with: Evacuators

#16660	White (50)
#16661	Blue (50)
#16662	Green (50)



VISUAL SYSTEMS & MISCELLANEOUS BOOK-KEEPING FORMS

#1662727	No. 1000 - Appt. & Daily Record	#1662878	No. 1 - Daily Efficiency Sheet
#1662731	No. 2000 - Exam Ledger Estimate	#1662883	No. 2B - Examination Sheet
#1662732	No. 2100 - Patient Ledger Only	#1662886	No. 3B - Service & Account Record
#1662737	No. 3600 - Disbursement & Exp Record	#1662899	No. 6000 - Patient Record Envelopes Buff (11.75 in. x 8.75 in.)
#1662711	No. 5000 - Yearly Recapitulation	#1662900	No. 6000 - Patient Record Envelopes Grey (11.75 in. x 8.75 in.)
#1662759	No. ADR-51 - Appointment & Daily Record (Undated)		

WONDER WEDGES®

Wonder Wedges are interproximal wooden wedges for holding matrix bands when performing all multi-surface fillings. Wonder Wedges keep filling and restorative materials in the matrix and reduce the formation of overhang, making polishing easier and quicker. Wonder Wedges are anatomically conformed (rounded tip and concave sides) to readily shape to the natural contour(s) of the mesial and distal spaces, minimizing the trauma to the patient's gingiva and decreasing discomfort. The end of the wedge is square for a secure grip. Made from sycamore wood, which swells slightly in a wet environment, Wonder Wedges fit tightly and will not loosen during dental treatment. Wonder Wedges are easily removed and will not chip or splinter. Choose from natural or a variety of colours in small, medium, large, or extra-large. Coloured wedges will not stain.

Tips: Select the largest wedge that will fit the interproximal space to prevent the wedge from slipping.

#0921417	12 mm Small, Natural (100)	#0921422	12 mm Small, Natural (500)
#0921418	13 mm Medium, Natural (100)	#0921423	13 mm Medium, Natural (500)
#0921419	15 mm Large, Natural (100)	#0921424	15 mm Large, Natural (500)
#0921420	17 mm Extra Large, Natural (100)	#0921425	17 mm Extra Large, Natural (500)
#0921421	Assorted, Natural (100)	#0921426	Assorted, Natural (500)

#0921480	12 mm Small, Blue (100)	#0921485	12 mm Small, Blue (500)
#0921481	13 mm Medium, Green (100)	#0921486	13 mm Medium, Green (500)
#0921482	15 mm Large, Yellow (100)	#0921487	15 mm Large, Yellow (500)
#0921483	17 mm Extra Large, Bleached (100)	#0921488	17 mm Extra Large, Bleached (500)
#0921484	Assorted, Colours (100)	#0921489	Assorted, Colours (500)



AMALGAM SEPARATION

- ☐ **Apavia® AVT Max**
Advanced Vacuum Technology
- ☐ **Apavia® HydroMax**
Custom Filtration System
- ☐ **Apavia® UltraMax**
POE Membrane Filtration System
- ☐ **Apavia® EcoMax**
Tap Water to Stabilized Ozone
- ☐ **Apavia® PurMax**
Drinking Water System



Apavia® is a fresh perspective in filtration technology. The Apavia® Drinking Water Systems are designed for simple and easy installation to improve the quality and taste of office water while being cost conscious and eco-friendly with a simple exchanging of the internal filter. The reusable Apavia® Drinking Water System housings minimize waste and maximize savings. Once the internal, long-life cartridge reaches maximum usage, the quarter turn release disconnect shuts off the water supply for easy maintenance and replacement. No spills and no lost time!

Apavia's AMB filters, made of a matrix of activated carbons, crystalline mono and di-calcium phosphates and chemically active calcium carbonate, all structurally bound into an interchangeable anti-microbial porous block. In other words, the filters absorb organic substances, heavy metals and pollutants to make your water as clean and great-tasting as possible!



Apavia® AVTMax

With the industry's highest retention rate of 99.7%, Apavia® provides an Advanced Vacuum Technology (AVT) solution that performs Amalgam Separation and other contaminant removal, protecting and extending the life of a vacuum pump. It not only captures mercury and other harmful particles, it protects and extends the life of the vacuum pump and drastically reduces the frequency of changing screen mesh traps. The Apavia® Amalgam Separator boasts a high containment capacity for extended service and pump life, cutting-edge advanced air circulation techniques with rapid velocity and waste capture. The AVTMax's fully moulded design offers a clear view of the system in action. With ISO 11143:2008 certification by NSF, meeting Best Management Practices (BMPs), and over 99.7% removal, Apavia's® high efficiency cyclonic filtration design leads the way in advanced vacuum solutions while also complying with the rigid amalgam separation requirement and all newly proposed EPA regulations.

The full AVTMax™ unit includes mounting bracket, AVTMax™ replacement canister and all installation parts to make the install a breeze with virtually no turnover time! Apavia's® large replacement canister is universally moulded to fit currently installed competitive units, allowing them to meet the requirements of the newly proposed regulations, while also extending the frequency of necessary change outs and reducing expenses. All replacement canisters are part of Apavia's® environmentally safe recycling programme that meets both state and federal regulations.

#39AVT



Apavia® EcoMax

Apavia® EcoMax™ technology turns ordinary tap water into stabilized ozone - a powerful natural cleaner, disinfectant, stain remover and deodourizer. After rapidly oxidizing all bacteria it comes in contact with in the 24-hour life cycle, EcoMax™ safely converts back into oxygen and water. The National Institute of Health has recorded that up to 65% of all microbial infections are caused by biofilm. Biofilm creates a shield protecting bacteria, cysts, fungi and parasites against most common disinfectants. The EcoMax™ system offers a powerful, naturally safe alternative that destroys the risk of these infection and disease-causing threats as well as disinfectant-resistant strains.

EcoMax™ technology turns ordinary tap water into stabilized aqueous ozone. It is safe to use on any surface, while eliminating the need for other cleaners, stain removers and deodourizers. With on-the-spot, continuous flow for spray bottles, mop buckets, carpet extractors and auto scrubbers, EcoMax™ kills germs and bacteria instantly on contact. EcoMax™ utilizes the oxygen that we breathe to create the most efficient, all-purpose cleaner available with 4,500 volts of electricity transforming O2 to O3. This combination creates a highly effective cleaning agent that breaks down dirt, grease, Bio-film and other contaminants in the same way as toxic cleaners, but naturally!

#39EMX

Apavia® HydroMax

The Apavia® HydroMax™ filtration system offers complete water filtration for your entire dental office. With multiple configuration options, the system is tailored for your office's specific water treatment needs, whether municipal or well water sources. If your water source has sulphur, irons, impurities, hardness, chlorine or just doesn't taste or smell right, HydroMax™ is the solution.

Determine your custom filter option to remove: Sediment, Chlorine, Chloramine, Pharmaceuticals, Taste & Odour, Sub-micron 0.2 filtration for cysts, heavy metals, Endotoxin, etc. In addition to providing your office with pure, filtered water, the HydroMax™ will help extend the life of office appliances (dishwashers, autoclaves, sanitizers, faucets, etc.), putting money back into your pocket!

Each HydroMax™ replacement filter comes with a unique handle designed top cap for lightweight and easy removal, a bag for proper disposal, and a double o-ring bottom connection into the housing for 1¼" PVC glue socket-plumbing connections. The patented no tool necessary housing allows for easy removal access of the filter for quick replacement, keeping your office running with no downtime.

#39HDM



Apavia® PurMax

The Apavia® Drinking Water Systems are designed for simple and easy installation to improve the quality and taste of your office water while being cost conscious and eco-friendly with a simple exchanging of the internal filter. The reusable Apavia® Drinking Water System housings minimize waste and maximize savings. Once the internal, long-life cartridge reaches maximum usage, the quarter turn release disconnect shuts off the water supply for easy maintenance and replacement.

The PurMax® systems connect directly to the cold water faucet for full flow filtration; providing your office instant, refreshingly clean water that is better for your equipment, dental lines and drinking needs. So ditch the bottles and jugs, and allow the PurMax® system to treat municipal and well water lines. Our AMB filters, are made of a matrix of activated carbons, crystalline mono and di calcium phosphates and chemically active calcium carbonate, all structurally bound into an interchangeable anti-microbial porous block. In other words, our filters adsorb organic substances, heavy metals and pollutants to make your water as clean and great-tasting as possible!

The PurMax® is ideal for the improvement of the unpleasant taste of chlorine and chloromine, while attacking metals, ammonia and other odors in the water lines of your dental water system. Activated carbon is recognized by the EPA as the best available control technology for many contaminants including: lead, chlorine, chloramine, pharmaceuticals, pesticides, arsenic and many others.

#39PMX



Apavia® UltraMax

The Apavia® UltraMax™ Filtration system uses a premiere barrier technology to protect your office from bacteria, virus, spores and cysts at the Point-of-entry (POE). The Apavia® UltraMax™ uses a membrane filtration process where water is pushed through a physical filter with a transmembrane pressure supplied by a pump or gravity. Water passes through the membrane pores and disperses effectively retaining suspended solids, viruses and bacteria on the outside of the membrane, preventing them from entering the service lines. The system is simple and safe to install, allowing it to physically prevent water-related diseases caused from microorganisms. It is well-known that treated municipal water systems and water derived from wells include bacteria and other potential health threats.

#39UMX



ACID ETCHANTS

Acid Etch

35% Acid Etch



WATCH THE
TECHNIQUE VIDEO!

ACID ETCH

Acid Etch is silica and glycerin-free, providing a cleaner, rougher, increased surface area free of bio-film. Acid Etch creates a deeper etch pattern for an easier application of bonding materials. For increased patient comfort and easy clean-up, Acid Etch is highly viscous, allowing for more accurate placement. Acid Etch is odourless, tasteless, and available in a 35% concentration aqua blue formula.

Bulk packaging available. Recommended for use with: Hydrox®

Tips: Acid Etch contains phosphoric acid which may cause burns or skin sensitization in susceptible individuals. Before using this product cover the patient's eyes with protective eyewear. Avoid contact with oral tissues, eyes and skin. Do not use Acid Etch on dental pulp tissue; cover dentin close to the pulp (< 1 mm remaining) with a hard setting calcium hydroxide liner (Bosworth Hydrox).

Gel d'acide phosphorique. Pour mordançage par acide de l'émail et la dentinePhosphorsäure-Ätzel. Zur Dentin- und Schmelzätzung • Gel di acido fosforico. Per mordenatura di smalto e dentinaGel de ácido fosfórico. Para cauterização a ácido de esmalte e dentina • Gel de ácido fosfórico. Para grabado ácido del esmalte y dentina • Fostorsyregel. För emalj- och dentinsyretsning

“ The Bosworth Acid Etch allows for flowability and clarity while providing complete coverage of the prepared surface. My patients prefer this easy flowing etchant compared to other etching gels due to this key feature. ”
Dr. Richard Bona

#0921346	Standard Kit	#0921345	1.2 ml Kit	#0921303	Jumbo Syringe Kit
4 - 3ml syringes Acid Etch 35% Phosphoric Acid Syringes		4 - 1.2ml syringes Acid Etch 35% Phosphoric Acid Syringes		50 ml syringe Acid Etch 35% plunger 3 - 3ml syringes (empty)	

ANAESTHETICS

- ComfortCaine™**
Topical Anaesthetic
- Gelato**
Topical Anaesthetic
- Gelato Benzocaine Anaesthetic Metered Spray**
Topical Anaesthetic
- Cora-Caine®**
Denture Analgesic Adhesive



COMFORTCAINE™

ComfortCaine is a topical anaesthetic which numbs the gums quickly for a short period (12-15 minutes) and contains 20% Benzocaine with virtually no systemic absorption. It is available as a gel in 1 oz wide mouth jars. ComfortCaine gel's smooth consistency lubricates as well as anaesthetizes.

Tips: ComfortCaine can also aid in the taking of impressions and intraoral radiographs for patients with excessive gag reflex.

In my dental practice, patient comfort is a top priority, which means using topical anaesthetic prior to administration of local anesthesia. ComfortCaine tastes great, comes in a variety of tasty flavours (my favourite is Pina Colada!) and helps make giving injections easier on the patient. I don't start a dental procedure without ComfortCaine. Dr. Robert A. Lowe

#16602	1 oz (30 g) Gel Piña Colada	#16605	1 oz (30 g) Gel Banana
#16603	1 oz (30 g) Gel Cherry	#16606	1 oz (30 g) Gel Grape
#16604	1 oz (30 g) Gel Strawberry	#16607	1 oz (30 g) Gel Mint



CORA-CAINE®

Cora-Caine is a pain-relieving adhesive ointment for new or intermediate dentures. The active ingredient is 16% Benzocaine. Cora-Caine provides immediate relief and a secure fit by soothing tender ridges and relieving the discomfort of sore spots. Acting as a gel-like cushion between the denture and the tissue, Cora-Caine is applied directly onto the denture to provide positive retention and enhance stability. Cora-Caine is pink in colour and has no smell or aftertaste. Available for the dental office in a box of 36 count 4 g tubes or for home use in a 1 oz tube. Remains effective up to 24 hours. Perfect product for denture break-in period.

Tips: Cora-Caine is a great product to include in a new denture kit. Do not use Cora-Caine for patients who have a history of allergy to local anaesthetics such as procaine, butacaine, or other “-caine” anaesthetics.

Dolor-Aliviando de dentadura Ungüento Adhesivo • Douleur-Soulager de dentier Baume Adhésif

#16621	36 - 4 g Tubes	case: 50	#16624	1 oz (28 g) Tube	case: 480
--------	----------------	----------	--------	------------------	-----------

GELATO TOPICAL ANAESTHETIC GEL

The Gelato 20% Benzocaine Anaesthetic Gels are fast-acting with no systemic absorption. The 20% benzocaine provides temporary relief of pain during procedures, including local anaesthetic injections, periodontal curettage, impression taking, scaling, intra-oral radiographs, root planning and prophylaxis. The gel comforts patients with ulcers, wounds or other minor mouth irritations. Gelato has no bitter aftertaste.

#03-02419	1 oz (.03 L) Bubble Gum	#03-01319	1 oz (.03 L) Grape
#03-02219	1 oz (.03 L) Strawberry	#03-02519	1 oz (.03 L) Piña Colada
#03-02319	1 oz (.03 L) Cherry	#03-02119	1 oz (.03 L) Mango
#03-03519	1 oz (.03 L) Banana	#03-02819	1 oz (.03 L) Mint
#03-02019	1 oz (.03 L) Raspberry	#03-01719	1 oz (.03 L) Piña Colada “Dye Free”



GELATO BENZOCAINE ANAESTHETIC METERED SPRAY

The Gelato Benzocaine Anaesthetic Spray is great for use on a large area that requires anaesthesia. Containing 20% benzocaine the spray works quickly and effectively – approximately 15 to 20 seconds.

#03-01519	Metered Spray (includes 25 dispensing tips)
#03-01619	Replacement Tips



DENTURE BASE, RELINE & REPAIR

BASE MATERIALS

- Impact 1500™**
Denture Resin Base
- Impact 2000™**
High-impact Denture Resin Base

DENTURE ANALGESIC ADHESIVE

- Cora-Caine®**
Denture Analgesic Adhesive

PRESSURE INDICATOR PASTE

- Hi-Spot®**
Pressure Indicator Paste

RELINE MATERIALS

- DentuSil™**
Silicone Relining Material (SOFT)
- Light Liner®**
Light-cure Chairside Reline Composite (LIGHT-CURE, HARD & SOFT)
- New Truliner®**
PEMA Denture Corrective Relining Material (HARD)
- Original Truliner®**
PMMA Denture Corrective Relining Material (HARD)
- Trusoft™**
Resilient Acrylic Denture Relining Material (SOFT)

REPAIR MATERIALS

- Duz-All®**
Acrylic Repair Material
- TruRepair®**
Fast-bonding Acrylic Repair Material
- Ultra Fast**
Repair Liquid

SURFACE POLISHES

- Glaze™**
Self-cure Resin Polish

TISSUE CONDITIONERS

- Softone®**
Denture Tissue Treatment


BASE MATERIALS

IMPACT 1500™

Impact 1500 is a denture base resin which provides stronger, more natural looking dentures with the added benefit of a 20-minute process time. Similar to Impact 2000, Impact 1500 provides the advantages of heat cure polymerization and self-curing. Impact 1500 will not warp or bleach during processing and finishing.

Bulk packaging available.

#092127	Standard Kit (Vein Powder)
---------	----------------------------



COMPONENTS	
#092128	1 lb (454 g) Powder Only (vein)
#092129	8 oz (236 ml) Liquid Only
#0921509	Impact Denture Bags (50)











H Ships Hazardous

IMPACT 2000™

Impact 2000 is a heat-cured denture base acrylic with high-impact strength and durability, formulated for 90 minute processing that will resist breakage and warping. Trims easily for chairside adjustments, Impact 2000 is batch controlled for colour stability and consistency. Impact 2000 powders are cadmium-free.

Bulk packaging available.

#0921500	Standard Kit
Specify powder shade: Vein, Dark Pink, or MP (Purple Pink)	
	
1lb. (454g) Powder	8 oz. (236ml) Liquid
	
20 - Denture Bags	
COMPONENTS	
#0921501	1 lb (454 g) Powder (Specify Shade)
#0921502	8 oz (236 ml) Liquid
#0921509	Impact Denture Bags (50)

#0921503	Bulk Kit
Specify powder shade: Vein, Dark Pink, or MP (Purple Pink)	
	
5lb. (2.3kg) Powder	1 qt. (946ml) Liquid
	
25 - Denture Bags	
COMPONENTS	
#0921504	5 lb (2.3 kg) Powder (Specify Shade)
#0921505	1 qt (32 oz, 946 ml) Liquid
#0921509	Impact Denture Bags (50)

DENTURE ANALGESIC ADHESIVE



CORA-CAINE®

Cora-Caine is a pain-relieving adhesive ointment for new or intermediate dentures. The active ingredient is 16% Benzocaine. Cora-Caine provides immediate relief and a secure fit by soothing tender ridges and relieving the discomfort of sore spots. Acting as a gel-like cushion between the denture and the tissue, Cora-Caine is applied directly onto the denture to provide positive retention and enhance stability. Cora-Caine is pink in colour and has no smell or aftertaste. Remains effective up to 24 hours. Perfect product for denture break-in period.

Tips: Cora-Caine is a great product to include in a new denture kit. Do not use Cora-Caine for patients who have a history of allergy to local anaesthetics such as procaine, butacaine, or other “-caine” anaesthetics.

Dolor-Aliviando de dentadura Ungüento Adhesivo • Douleur-Soulager de dentier Baume Adhésif

#16621	36 - 4 g Tubes	case: 50	#16624	1 oz (28 g) Tube	case: 480
--------	----------------	----------	--------	------------------	-----------

PRESSURE INDICATOR PASTE



HI-SPOT®

Hi-Spot is a quick and economical method for determining sore spots and high-pressure areas. Apply to dentures, bridges and inlay cavities for inspection against binding. Hi-Spot's smooth and creamy texture makes it easy to work with and it is tasteless and odourless. Great for implant adjustments.

#0921430	Standard Kit	#0921440	Bulk Kit
			
2 - 14g Paste		4 - 14g Paste	

RELINE MATERIALS

DENTUSIL™

DentuSil is a cartridge dispensed silicone soft reline material that provides denture patients with a comfortable, reliable reline in one office visit. Flexible and resilient, DentuSil's silicone base has a low porosity, making it resistant to stains and odours. The translucent pink colour makes for a natural, healthy appearance. Adhesive holds the material in place, securing the reline in minutes. Contains approximately 10 applications. A successful reline with proper maintenance will last one year. Clean with a soft brush.

Tips: DentuSil contains no monomer and is perfect for use on patients with sensitivity to acrylics. DentuSil can be used to cushion dentures as overdentures heal. Use DentuSil in the dental laboratory by placing the denture on the patient's model in the dental flask. Do not use with condensation reaction type silicone materials or adhesive.

Recommended for use with: E-Z Squeeze™ Universal Dispenser and Static Mixing Tips for 5 inch cartridges.




COMPONENTS	
#0921276	Standard Kit
	
50 ml Cartridge	Mixing Tips
	
10 ml Liquid Adhesive	10 ml Silicone Release

LIGHT LINER®

Light Liner is a chairside reline material that light cures in minutes. Available in a soft or hard non-methyl methacrylate formula, this microfilled composite flows easily and offers the convenience of a one visit procedure without multiple pour-ups or impressions. Accurate relines and quality repairs are achieved in minutes at a fraction of the cost of sending it to the lab. Light Liner is a dual-cure, self-bonding material that will not shrink, resists staining, and is colour stable. Light Liner is virtually as strong as the denture, chemically bonding to it to prevent flaking or peeling. FOR DENTURE RELINING ONLY.

Recommended for use with: Cora-Caine® and Masque™ | Bulk packaging available.

Composite de regarnissage dentaire photopolymérisable • Lichthärtender Unterfütterungskomposit • Composito per ribasatura polimerizzabile per applicazioni Chairside • Composto de revestimento imediato fotopolimerizável • Compuesto de revestimiento Chairside para curar con luz

COMPONENTS	
#0921463	Standard Kit, Hard
#0921464	Standard Kit, Soft
Packaged in 8 single-dose packs for convenience.	
	
8 - 6g Composite Powder	8 - 4g Composite Liquid
	
	30cc Air Barrier

NEW TRULINER®

New Truliner is a PEMA denture corrective relining material. New Truliner's advanced formula reduces stinging and burning sensation which allows for a restoration on a denture's original fit in as little as 5 minutes. Colour stability and a high density surface guarantee greater patient satisfaction. New Truliner chemically fuses to the denture in 15-20 minutes and is formulated to last up to a year. Once cured, rebased surfaces can be trimmed and polished like the original denture without showing a line of demarcation.




Mixing Ratio: 7.6 ml liquid (1 full glass vial) to 15 ml powder (1 full plastic vial).

Tips: Before making a permanent denture reline, treat the tissue with a tissue conditioner such as Softone to eliminate the possibility of the final corrective reline duplicating the disturbed or deformed mucosa caused by the deficient denture.

Recommended for use with: Cora-Caine®, Glaze™, Gorilla Seal™, Masque™, and Pressure Pot | Bulk packaging available.

Correction de prothèse matériau de remodelage • Gebißkorrektur-unterfütterungsmaterial • Materiale di rivestimento per la correzione di dentiere • Material de revestimento corretivo para prótese • Material de revestimento correctivo de prótesis dental

#0921970	Standard Kit, Pink	case: 14
#0921971	Standard Kit, Vein	case: 14
#0921972	Standard Kit, Clear	case: 14

		
8 oz. (227g) Powder	6 oz. (177ml) Liquid	1/3oz. (10ml) Bonding Liquid

COMPONENTS	
#0921973	8 oz (227 g), Powder, Pink
#0921974	8 oz (227 g) Powder, Vein
#0921975	8 oz (227 g) Powder, Clear
#0921976	6 oz (177 ml) Liquid
#0921983	1/3 oz (10 ml) Bonding Liquid

#0921977	1 lb (454 g) Bulk Kit, Pink
#0921978	1 lb (454 g) Bulk Kit, Vein
#0921979	1 lb (454 g) Bulk Kit, Clear

		
1lb. (454g) Powder	8oz. (236ml) Liquid	1/3oz. (10ml) Bonding Liquid

COMPONENTS	
#0921980	1 lb (454 g) Powder, Pink
#0921981	1 lb (454 g) Powder, Vein
#0921982	1 lb (454 g) Powder, Clear
#0921986	8 oz (236 ml) Liquid
#0921984	4 oz (118 ml) Bonding Liquid



H Ships Hazardous

Note: New Truliner liquid and powder is NOT interchangeable with Original Truliner liquid and powder. New Truliner is a PEMA formula and Original Truliner is a PMMA formula.



Note: New Truliner liquid and powder is NOT interchangeable with Original Truliner liquid and powder. New Truliner is a PEMA formula and Original Truliner is a PMMA formula.

ORIGINAL TRULINER®

Original Truliner is a PMMA colour stable, permanent reline material formulated for extended life (3-5 years). The extremely fine-milled powder cures to a dense consistency and prevents the reline from assuming odours, tastes, or stains. Original Truliner is a unique reverse curing acrylic (cures from the inside out) which releases the free monomer before it is inserted in the mouth, eliminating chemical and thermal sensitivity and forming a permanent bond in as little as 10 minutes that will not chip or discolour. Can also be used for repairs and orthodontic appliances.


Mixing Ratio: 7.6 ml liquid (1 full glass vial) to 15 ml powder (1 full plastic vial).

Tips: Before making a permanent denture reline, treat the tissue with a tissue conditioner such as Bosworth Softone to eliminate the possibility of the final corrective reline duplicating the disturbed or deformed mucosa caused by the deficient denture.

Recommended for use with: Cora-Caine®, Glaze™, Gorilla Seal™, Masque™, and Pressure Pot | Bulk packaging available.

Correction de prothèse matériau de remodelage • Gebilkkorrektur-unterfütterungsmaterial • Materiale di rivestimento per la correzione di dentiere • Material de revestimento correctivo para prótese • Material de revestimiento correctivo de prótesis dental


#0921940	Standard Kit, Pink	case: 20
#0921941	Standard Kit, Vein	case: 20
#0921942	Standard Kit, Clear	case: 20



COMPONENTS

#0921943	8 oz (227 g) Powder, Pink
#0921944	8 oz (227 g) Powder, Vein
#0921945	8 oz (227 g) Powder, Clear
#0921946	6 oz (177 ml) Liquid

#0921951	1 lb (454 g) Bulk Kit, Pink
#0921952	1 lb (454 g) Bulk Kit, Vein
#0921953	1 lb (454 g) Bulk Kit, Clear



COMPONENTS

#0921947	1 lb (454 g) Powder, Pink
#0921948	1 lb (454 g) Powder, Vein
#0921949	1 lb (454 g) Powder, Clear
#0921950	8 oz (236 ml) Liquid




TRUSOFT™

Trusoft is a plasticized acrylic PEMA that is specifically designed as an intermediate liner for patients requiring a prolonged adjustment period prior to using a hard liner or denture base material (Softone is for shorter adjustment periods). Its elasticity will ensure against pressure pain, tender ridges, sensitive areas, and irritated tuberosities. Trusoft is formulated to remain resilient for up to 30 days. Can be used for every patient as a precautionary measure but is especially necessary for patients who are slow to heal or sensitive to wearing a new denture. Patient is ready for a permanent reline (New or Original Truliner) when he or she can exhibit control and acceptance with Trusoft.

Mixing Ratio: 6.8 ml (glass vial) to 9 g (plastic vial) powder. | Tips: Trusoft can be used in between tissue treatments with Softone. Recommended for use with: CoraCaine® and Masque™ | Bulk packaging available.

Matériau de remodelage élastique • Elastisches Unterfütterungsmateria • Materiale per rivestimento resiliente • Materiale acrílico elástico per ribasatura intermedia • Resina Suave para Rebasado Intermedio de Dentaduras


#0921250	Standard Kit
----------	--------------



COMPONENTS

#0921251	3 oz (85 g) Powder
#0921252	4 oz (118 ml) Liquid
#0921983	1/3 oz (10 ml) Bonding Liquid

#0921253	Bulk Kit
----------	----------



COMPONENTS & BULK

#0921254	1 lb (454 g) Powder
#0921255	8 oz (236 ml) Liquid
#0921983	1/3 oz (10 ml) Bonding Liquid
#0921984	4 oz (118 ml) Bonding Liquid

REPAIR MATERIALS




DUZ-ALL®

A PMMA all-purpose, self-curing acrylic repair material intended for any emergency including: replacing broken teeth, repairing fractures, adding clasps, making individual impression trays, cold-cure full dentures, orthodontic retainers, rebasing and correcting full or partial dentures, and adding a post-dam and appliances. Sets in less than 20 minutes with no special equipment needed. Ideal for making or repairing jacket crowns.

Bulk packaging available.

#166260	Standard Kit
Specify Pink, White, Clear, or Vein powder.	



COMPONENTS & BULK

#166261	8 oz (227 g) Powder (Specify Shade)
#166262	4 oz (118 ml) Liquid
#166264	1 lb (454 g) Powder (Specify Shade)
#166265	8 oz (236 ml) Liquid Only
#166269	32 oz (946 ml) Liquid

TRUREPAIR®


TruRepair is a fast-bonding acrylic repair material for repairing fractured dentures, replacing broken or dislodged denture teeth, and any additional office or lab emergency. A special internally-pigmented spherical globe of powder blends with denture base materials to produce a colour-stable match. TruRepair liquid is cross-linked.

Mixing Ratio: 2-3 parts powder to 1 part liquid.

Recommended for use with: Pressure Pot and Masque™ | Bulk packaging available.

#092160	Standard Kit
---------	--------------

Specify Pink, Vein, or Clear Powder.



COMPONENTS

#092161	3 oz (85 g) Powder, Pink
#092162	3 oz (85 g) Powder, Vein
#092163	3 oz (85 g) Powder, Clear
#092164	4 oz (118 ml) Liquid

ULTRA FAST

Ultra Fast is a cross-linked repair liquid that can be used with any acrylic powder to produce a self-curing denture material. Use Ultra Fast with any self-curing acrylic liquid to accelerate setting time.

#0921355	8 oz (236 ml) Bottle
----------	----------------------

#0921357	32 oz (946 ml) Can
----------	--------------------

SURFACE POLISH

GLAZE™

Glaze is a fast-setting, self-cure resin which reduces polishing time and dries to a clear, hard surface with a high gloss finish. Glaze can be brushed onto a clean composite or acrylic material to achieve a polished surface that enhances the appliance and reduces the need for additional polishing. Recommended for use on provisional, reline, and orthodontic appliances. Glaze fills microscopic voids and reduces bacterial growth. Sets in 20 seconds.

Finierlack Für die Orthodontie, Zahnprothesen, Kronen, Unterfütterungen und Provisorien • Glaçure de finissage Pour appareils d'orthodontie, prothèses, couronnes, rebasages et restaurations temporaires • Smalto per finitura Per apparecchi ortodontici, dentiere, capsule, ribasature e provvisori • Verniz de acabamento Para aparelhos ortodônticos, dentaduras, coroas, revestimentos e provisórias • Esmalte para acabado Para aparatos ortodônticos, prótesis dentales, coronas, revestimientos y aplicaciones temporales

#0921236	Standard Kit
----------	--------------



TISSUE CONDITIONERS

SOFTONE™

Softone is a resilient acrylic denture tissue conditioner and functional impression material that conditions, massages, and rehabilitates damaged tissue. Softone's unique flow and natural give stimulate the natural healing process and restores tissue to a sound physiological state. Softone has a clean, refreshing taste that will not absorb odours for greater patient comfort and compliance. Softone allows the denture to stay in the mouth throughout the treatment process. For patients requiring a temporary soft reline or prolonged adjustment period before the permanent hard reline, use Bosworth Trusoft.

Mixing Ratio = 6.8 ml (1 glass vial) to 9.0 g (1 plastic vial) of powder.


Tips: Polishing may be expedited by first chilling denture in refrigerator or ice water. Softone is a temporary impression material formulated to last up to 4 days in the mouth.

Recommended for use with: TAC™ | Bulk packaging available.

Traitement pour tissu de prothèse dentaire souple à l'acrylique • Gewebebehandlung bei dehnbaren Zahnprothesen aus Acryl • Trattamento tessuto per protesi in acrilico resiliente • Tratamento dos tecidos, prótese dentária acrílica maleável • Tratamiento acrílico resistente para tejido de dentaduras

#0921775	Standard Kit	case: 20
----------	--------------	----------

Specify Pink or White Powder.




COMPONENTS

#0921776	3 oz (85 g) Powder (Specify Pink or White)
#0921777	4 oz (118 ml) Liquid

#0921778	Bulk Kit
----------	----------

Specify Pink or White Powder.



COMPONENTS

#0921779	1 lb (454 g) Powder (Specify Pink or White)
#0921780	8 oz (236 ml) Liquid



HYGIENE

- ☐ **Simplastic® Dispenser**
Acrylic Dispensing Unit
- ☐ **Simplastic® Hazard Bag**
10-gallon Hazard Waste Disposal Bag
- ☐ **Simplastic® Patient Drape**
Full Coverage Patient Drape
- ☐ **Simplastic® Sleeve**
Barrier Sleeve 4 in x 6 in
- ☐ **Simplastic® Tray Cover**
Multi-Purpose Cover 12 in x 17 in
- ☐ **Simplastic® Unicover**
Multi-Purpose Gusseted Barrier
16.5 in x 13.5 in x 31 in
- ☐ **Vigilance**
Multi-purpose Enzyme Concentrate Solution



SIMPLASTIC® BARRIER PROTECTIVE SYSTEM

- SIMPLASTIC DISPENSER**
A clear, acrylic dispensing unit for all Simplastic products. Can be wall-mounted.
Dimensions: 11.5 in x 17 in x 10.5 in (29.2 cm x 43.2 cm x 26.7 cm)
- SIMPLASTIC HAZARD BAG**
A 10-gallon (37.9 L) waste bag specifically designed for the disposal of infectious and hazardous waste.
- SIMPLASTIC PATIENT DRAPE**
Full drape which provides total coverage for patients from debris. Controls cross contamination from patients' skin and clothing.
Dimensions: 45 in x 48 in (114.3 cm x 122 cm)
- SIMPLASTIC SLEEVE**
A barrier sleeve for protecting curing lights, three-way syringes, and evacuation tubes. Easily dispensed from the uni-bag handles. Each handle holds 100 sleeves.
Dimensions: 4 in x 16 in (10.2 cm x 40.6 cm)
- SIMPLASTIC TRAY COVER**
A multi-purpose cover for protecting trays, instrument carts, keyboards, and film. Easily dispensed from the uni-bag handles. Each handle holds 100 covers.
Dimensions: 12 in x 17 in (30.5cm x 43.2cm)
- SIMPLASTIC UNICOVER**
An anti-static gusseted barrier for protecting chairs, lights, and x-ray heads. Easily dispensed from rolls or from the Simplastic Dispenser.
Dimensions: 16.5 in x 13.5 in x 31in (41.9 cm x 34.3 cm x 78.7 cm)

#0921673	Dispenser, Clear Acrylic	#0921676	Hazard Bag (100)
#0921662	Unicover (300 Roll)	#0921660	Sleeve (1000)
#0921670	Unicover (650 Roll)	#0921671	Sleeve (2000)
#0921661	Tray Cover (1000)	#0921674	Patient Drape (400)
#0921672	Tray Cover (2000)		

VIGILANCE®

Vigilance is a proteolytic enzyme detergent specially formulated for removing bioburden and contaminants from instruments, working surfaces, sink traps and clothing. Its pleasing aroma, neutral pH, and reduced foaming action make it the ideal ultrasonic solution, evacuation cleaner, and instrument pre-soak. Vigilance meets hospital standards for a detergent wash, contains a rust inhibitor, and is biodegradable, non-corrosive, and environmentally safe.

Dilution Mixing Ratio: 1 to 2 ounces of Vigilance per gallon of water. | Tips: Shake well before use.

#0921469	1/2 Gallon (1.9 L) - Calibrated	#0921472	1 Gallon (3.8 L) Refill	case: 4
#0921474	Pump attachment	#0921474	Pump attachment	



IMPRESSION

- ADHESIVE**

☐ **TAC™**
Tray Adhesive Compound

☐ **Big Bite Trays™**
Disposable Bite Registration Trays

☐ **Re'CORD™**
Polyvinylsiloxane Bite Registration Material

☐ **Superbite™**
Zinc Oxide Eugenol Bite Registration Paste
- DISPOSABLE & STAINLESS STEEL IMPRESSION TRAYS**

☐ **BoTrays®**
Disposable Impression Trays

☐ **Edentulous Tray Aways®**
Edentulous Disposable Impression Trays

☐ **Stainless Steel Impression Trays**
100% Stainless Steel Impression Trays

☐ **Tray Aways®**
Disposable Impression Trays
- IMPRESSION MATERIALS**

☐ **Plastopaste™**
Zinc Oxide Eugenol Impression Paste (Light Body)

☐ **Rimseal™**
Acrylic Impression Material (Yellow)

☐ **Sapphire™**
Acrylic Impression Material (Blue, Clear)

☐ **Supergel®**
Alginate Impression Material

☐ **Superpaste™**
Zinc Oxide Eugenol Impression Paste (Heavy Body)

☐ **Supersil™**
Polyvinylsiloxane Impression Material
- CLEANERS**

☐ **G.T.C.™**
General Tray Cleaner

☐ **Orange Solvent**
Wax and Impression Material Solvent
- CUSTOM TRAY MATERIALS**

☐ **Fastray™**
Custom Tray and Acrylic Base Plate Material

☐ **Fastray™ LC**
Light-cure Custom Tray and Acrylic Base Plate Material
- PLASTER**

☐ **Plastogum®**
Impression Corrective Plaster

ADHESIVE


TAC™

TAC is a tray adhesive compound available in a convenient spray can or brush-on liquid for use with alginate impression materials. TAC securely holds the alginate in the tray when the impression is removed from the patient's mouth. Excellent for orthodontic use.

Recommended for use with: Edentulous Tray Aways® and Tray Aways®.

Tips: Spray TAC on denture before applying soft reline to facilitate easy removal. Clean up with 70% alcohol.

Composé adhésif pour plateau pour l'alginate • Klebende Zusammensetzung für abformlöffeln für alginat • Adesivo compostoper vassoio per alginato • Composto adesivo para impressão em moldura • Compuesto adhesivo para cubetas para el alginato

#0921881	Standard Kit (Brush-on Liquid)	case 150
 2 - 0.5 oz (14ml) Bottles		
#0921882	4 oz (118 ml) Bulk Bottle	

#0921880	Spray - 3.5 oz (99 g)	case 60
----------	-----------------------	---------



H Ships Hazardous

BITE REGISTRATION




BIG BITE TRAYS™

Big Bite Trays are economical, disposable bite registration trays with plastic frames and a gauze insert for use when recording centric occlusion. Big Bite Trays can be used on either side to cover a quadrant or for bi-lateral registration by clipping together right and left trays. The thin gauze mesh does not interfere with closure and ensures an accurate bite.

Recommended for use with: Re'CORD™, Superbite™

Dimensions:
A width 28.5 mm
B height 44.54 mm
C depth 9.5 mm



#0921841 50 Trays (25 left, 25 right)




RE'CORD™

Re'CORD is a polyvinylsiloxane bite registration material specially formulated to provide enhanced accuracy, convenience, and patient acceptance. Re'CORD is blue in colour, odourless, and tasteless. A fast-setting material, Re'CORD achieves snap-set hardness in under 1 minute to reduce patient discomfort and prevent costly remakes due to patient reflexive movement. Directly apply Re'CORD to the occlusal surface or use with a Bosworth Big Bite Tray for effortless bite registration in centric occlusion or at rest. Re'CORD remains virtually distortion free weeks after use.

Recommended for use with: E-Z Squeeze™ Universal Dispenser and Big Bite Trays™

		COMPONENTS
#0921315	Standard Kit	#0921414 Static Mixing Tips for 5 inch cartridges (48)
#0921316	Standard Kit (Rigid)	



2 - 50ml Cartridge
10 - Mixing Tips




SUPERBITE™

Superbite is a traditional zinc oxide eugenol paste that combines controlled flowability and non-sagging body to enable unexcelled reproduction of the oral tissues that support a denture, registration of centric, and accurate seating of the casts for articulation. Superbite paste is a blue, two-component material, with a 1 minute in the mouth setting time. The absence of secondary flow and resistance to closure will enhance patient comfort.

Recommended for use with: Big Bite Trays™

Zinkoxyd-Eugenol-Bißabdruckpaste • Pasta para impresión de mordida a base de eugenol de óxido de zinc • Pâte pour prise d'empreintes à l'eugénol d'oxyde de zinc • Pasta gelificante all'eugenolo mescolato con ossido di zinco per laregistrazione delle impronte dentarie

		COMPONENTS
#0921815	Standard Kit	#0921816 Tube Base
		#0921817 Tube Catalyst



7.5oz. (212g) Base
3.25oz. (92g) Catalyst

CLEANERS

G.T.C.™

G.T.C. is a special oxygenating formula designed to clean and brighten trays (including stainless steel trays), instruments, oral evacuators, diamond instruments and burs while removing alginate, zinc oxide eugenol paste, wax, and modelling compounds. G.T.C. contains a special flavouring agent and is effective as a deodourizer. G.T.C. is excellent for removing dried blood and for flushing and cleaning evacuators.

Mixing Ratio: 2 oz (60 ml) of G.T.C. to 1 qt of warm water. | Bulk packaging available.

#0921459 2 lb (908 g) Bottle

case: 20

#0921004 100 lb (45.45 Kg) Drum

ORANGE SOLVENT

Orange Solvent is a wax and impression material solvent that quickly removes wax, baseplate, and resinous impression materials from acrylic teeth and denture bases. When used according to directions, Orange Solvent can be used on face and hands without injury to the most delicate skin. Orange Solvent does not contain carbon tetrachloride, acetone, chloroform, or any solvent which will injure, dissolve, or craze acrylic teeth or bases. Apply only with cotton or gauze.

Bulk packaging available.

#0921575 16 oz (474 ml)

CUSTOM TRAY MATERIALS


FASTRAY™

Fastray is a dimensionally stable, radiopaque, self-curing individual custom tray and acrylic base plate material. Fastray has no memory and when used with a separator, it will not stick to fingers or the model.

Recommended for use with: Masque™ | Mixing Ratio: 1 oz (26 g) powder to 6 ml of liquid. | Bulk packaging available.

Porte-empreinte individuel et matériau pour plaque-base acrylique • Autopolymerisat für individuelle Abformlöffel und Kunststoff-Basisplatten • Materiale per portaimpronte individuali e per basette in resina acrilica • Material para placas de base acrílicas e moldes individuais • Material para cubetas individuales y bases acrílicas


#0921375	Standard Kit, Regular Blue	case: 50
#0921376	Standard Kit, Regular Pink	case: 50
#0921377	Standard Kit, Regular White	case: 50
#0921378	Standard Kit, X-Fast Blue	case: 50
#0921379	Standard Kit, X-Fast Pink	case: 50
#0921380	Standard Kit, X-Fast White	case: 50



1lb. (454g) Powder
8 oz. (236ml) Liquid

COMPONENTS & BULK	
#0921381	1 lb (454 g) Powder Blue
#0921382	1 lb (454 g) Powder Pink
#0921383	1 lb (454 g) Powder White
#0921384	8 oz (236 ml) Liquid Regular
#0921385	8 oz (236 ml) Liquid X-Fast
#0921429	25 lb (11.34 kg) Powder (specify colour)

#0921386	Bulk Kit, Regular Blue
#0921387	Bulk Kit, Regular Pink
#0921388	Bulk Kit, Regular White
#0921389	Bulk Kit, X-Fast Blue
#0921390	Bulk Kit, X-Fast Pink
#0921391	Bulk Kit, X-Fast White



5lb. (2.2kg) Powder
32oz. (946ml) Liquid

COMPONENTS & BULK	
#0921392	5 lb (2.2 kg) Powder Blue
#0921393	5 lb (2.2 kg) Powder Pink
#0921394	5 lb (2.2 kg) Powder White
#0921395	1 qt. (32 oz, 946 ml) Liquid Regular
#0921396	1 qt. (32 oz, 946 ml) Liquid X-Fast
#0921415	1 gal (3.8 ml) Liquid Regular
#0921416	1 gal (3.8 ml) Liquid X-Fast


FASTRAY™ LC

Fastray LC is a light-cured custom tray material supplied in pre-formed sheets. When used with a separator, Fastray LC will not distort or stick to fingers and models. Procedure time is reduced and odour is eliminated due to no powder and liquid mixing. Fastray LC provides a firm, stable custom tray that is ready to use immediately after trimming and polishing.

Recommended for use with: Masque™ | Tips: Fastray LC is a light sensitive material. Keep box closed. Material may set after 15 minutes of light exposure.

Porte-empreinte individuel et plaque-base en acrylique • Lichtgehärtetes Autopolymerisat für individuelle Abformlöffel und Kunststoff-Basisplatten • Materiale fotopolimerizzabile per portaimpronte personalizzate e basi provvisorie per protesi acriliche • Material para placas de base acrílicas e moldes individuais fotopolimerizável • Material Fotocurable para Cubetas Individuales y Bases Acrílicas

#0921405	Standard Kit, Blue (50 Wafers)	case: 10
#0921406	Standard Kit, Pink (50 Wafers)	case: 10
#0921407	Standard Kit, White (50 Wafers)	case: 10



H Ships Hazardous



H Ships Hazardous



WATCH THE
TECHNIQUE VIDEO!

DISPOSABLE & STAINLESS STEEL IMPRESSION TRAYS



BOTRAYS®

BoTrays disposable impression trays are an economical choice for consistent, unsurpassed impressions. Alginate “grippers” and slotted walls keep material in the tray and reduce gag effect. Made from rigid, high-impact FDA approved styrene and organic FDA approved dye.

Recommended for use with: TAC, Supergel® and Supersil™

A = Width, B = Height, C = Depth

#1 Perforated Large Upper	A: 3.06", 77.6 mm	B: 2.41", 61.2 mm	C: 0.88", 22.4 mm
#2 Perforated Large Lower	A: 3.02", 76.7 mm	B: 2.30", 58.3 mm	C: 0.83", 21.1 mm
#3 Perforated Medium Upper	A: 2.91", 73.9 mm	B: 2.21", 56.2 mm	C: 0.83", 21.3 mm
#4 Perforated Medium Lower	A: 2.93", 74.3 mm	B: 2.18", 55.2 mm	C: 0.79", 20.0 mm
#5 Perforated Small Upper	A: 2.66", 67.6 mm	B: 2.00", 50.9 mm	C: 0.86", 21.9 mm
#6 Perforated Small Lower	A: 2.78", 70.6 mm	B: 2.06", 52.3 mm	C: 0.77", 19.6 mm
#7 Perforated Partial, UR, LL	A: 2.24", 56.9 mm	B: 2.29", 58.3 mm	C: 0.79", 20.1 mm
#8 Perforated Partial, UL, LR	A: 2.15", 54.5 mm	B: 2.33", 59.3 mm	C: 0.86", 21.8 mm
#9 Perforated Lower Anterior	A: 2.56", 65.0 mm	B: 1.50", 38.1 mm	C: 0.75", 19.0 mm
#10 Perforated Upper Anterior	A: 2.52", 64.1 mm	B: 1.47", 37.2 mm	C: 0.76", 19.2 mm

#0921364	#1 Perforated Large Upper (12)
#0921365	#2 Perforated Large Lower (12)
#0921366	#3 Perforated Medium Upper (12)
#0921367	#4 Perforated Medium Lower (12)
#0921368	#5 Perforated Small Upper (12)

#0921369	#6 Perforated Small Lower (12)
#0921370	#7 Perforated Partial Upper Right, Lower Left (12)
#0921371	#8 Perforated Partial Upper Left, Lower Right (12)
#0921372	#9 Perforated Lower Anterior (12)
#0921373	#10 Perforated Upper Anterior (12)



Clinicians Report®
August 2014
Evaluation

EDENTULOUS TRAY AWAYS®

Edentulous Tray Aways® are single-use, disposable perforated edentulous impression trays that can be easily drilled for implant impressions. These unique and economical trays allow for maximum material retention and accuracy without the use of an adhesive. These high impact styrene trays provide an excellent impression with outstanding details. Edentulous Tray Aways are anatomically designed to be patient friendly and to allow for a perfect fit in just one visit.

Recommended for use with: TAC™, Supergel®, and Supersil®

#0921600	#1 Upper Small (12)
#0921599	#2 Lower Small (12)
#0921700	#3 Upper Medium (12)
#0921701	#4 Lower Medium (12)

#0921702	#5 Upper Large (12)
#0921703	#6 Lower Large (12)
#0921706	#7 Upper X-Large (12)
#0921707	#8 Lower X-Large (12)

#0921704	Variety Bag (2 of each size)
----------	------------------------------



TRAY AWAYS®

Tray Aways are the disposable impression trays of choice suited for use with alginate or other types of impression materials where detailed reproduction and dimensionally-stable results are vital. Made from rigid, high-impact FDA approved styrene, Tray Aways assist in generating a consistent, unsurpassed impression free of distortion or bite through. Smooth, rounded edges, interior retention slots, and vertical perforations enhance material retention.

Recommended for use with: TAC™ and Supergel®.

A = Width, B = Height, C = Depth

#1 Perforated Large Upper	A: 3.06", 77.6 mm	B: 2.41", 61.2 mm	C: 0.88", 22.4 mm
#2 Perforated Large Lowe	A: 3.02", 76.7 mm	B: 2.30", 58.3 mm	C: 0.83", 21.1 mm
#3 Perforated Medium Upper	A: 2.91", 73.9 mm	B: 2.21", 56.2 mm	C: 0.83", 21.3 mm
#4 Perforated Medium Lower	A: 2.93", 74.3 mm	B: 2.18", 55.2 mm	C: 0.79", 20.0 mm
#5 Perforated Small Upper	A: 2.66", 67.6 mm	B: 2.00", 50.9 mm	C: 0.86", 21.9 mm
#6 Perforated Small Lower	A: 2.78", 70.6 mm	B: 2.06", 52.3mm	C: 0.77", 19.6 mm
#7 Perforated/Solid Partial, UR, LL	A: 2.24", 56.9 mm	B: 2.29", 58.3 mm	C: 0.79", 20.1 mm
#8 Perforated/Solid Partial, UL, LR	A: 2.15", 54.5 mm	B: 2.33", 59.3 mm	C: 0.86", 21.8 mm
#9 Perforated/Solid Lower Anterior	A: 2.56", 65.0 mm	B: 1.50", 38.1 mm	C: 0.75", 19.0 mm
#10 Perforated/Solid Upper Anterior	A: 2.52", 64.1 mm	B: 1.47", 37.2 mm	C: 0.76", 19.2 mm

Porte-Empreinte Jetable • Einweg-Abdruck-Tray • Vassio Per Impronte Monouso • Moldeira Descartável Para Impressão • Bandeja Desechable Para Impresiones

#0921884	#1 Perforated Large Upper (12)
#0921887	#2 Perforated Large Lower (12)
#0921885	#3 Perforated Medium Upper (12)
#0921888	#4 Perforated Medium Lower (12)
#0921886	#5 Perforated Small Upper (12)
#0921889	#6 Perforated Small Lower (12)
#0921890	#7 Solid Partial Upper R, Lower L (12)

#0921891	#8 Solid Partial Upper L, Lower R (12)
#0921892	#9 Solid Lower Anterior (12)
#0921893	#10 Solid Upper Anterior (12)
#0921894	#7s Perforated Partial Upper R, Lower L (12)
#0921895	#8s Perforated Partial Upper L, Lower R (12)
#0921896	#9s Perforated Lower Anterior (12)
#0921897	#10s Perforated Upper Anterior (12)

100% STAINLESS STEEL IMPRESSION TRAYS

Bosworth's Stainless Steel Impression Trays are created for use with alginate or any other type of impression material where detailed reproduction and dimensionally stable results are vital. They lock impression material in the tray to prevent distortion when removing it from the mouth. Uniform thickness gives superior rigidity for consistent results. Bosworth Stainless Steel Impression Trays have a smooth finish to minimize tissue trauma, and are thin enough to allow for optimum patient comfort. They are autoclavable for extended use and wear, and come in solid or perforated styles in a large assortment of pedodontic, edentulous, regular, and XL sizes for partial or full impressions.

Twist trays are the same great impression tray with a maneuverable handle to fit both upper and lower quadrants. Just spin the handle, and place in desired area - upper right, upper left, lower right, or lower left! Autoclavable for extended use and wear. Available in either solid or perforated.



#0921431	#1 Edentulous Solid Upper
#0921432	#2 Edentulous Solid Upper
#0921433	#3 Edentulous Solid Upper
#0921434	#4 Edentulous Solid Upper
#0921435	#5 Edentulous Solid Lower
#0921436	#6 Edentulous Solid Lower
#0921437	#7 Edentulous Solid Lower
#0921438	#8 Edentulous Solid Lower
#0921335	Set 104 (8) 4 Upper, 4 Lower

#0921337	#1 - Edentulous Perforated Upper
#0921338	#2 - Edentulous Perforated Upper
#0921339	#3 - Edentulous Perforated Upper
#0921340	#4 - Edentulous Perforated Upper
#0921341	#5 - Edentulous Perforated Lower
#0921342	#6 - Edentulous Perforated Lower
#0921343	#7 - Edentulous Perforated Lower
#0921344	#8 - Edentulous Perforated Lower
#0921307	Set 103 (8) - 4 Upper, 4 Lower

#0921441	#1 Regular Solid Upper
#0921442	#2 Regular Solid Upper
#0921492	#3 Regular Solid Upper
#0921493	#4 Regular Solid Upper
#0921443	#5 Regular Solid Upper
#0921444	#1 Regular Solid Lower
#0921445	#2 Regular Solid Lower
#0921494	#3 Regular Solid Lower
#0921495	#4 Regular Solid Lower
#0921446	#5 Regular Solid Lower
#0921440	Set 105 (10) 5 Upper, 5 Lower

#0921654	#1 Regular Perforated Upper
#0921668	#2 Regular Perforated Upper
#0921475	#3 Regular Perforated Upper
#0921476	#4 Regular Perforated Upper
#0921663	#5 Regular Perforated Upper
#0921664	#1 Regular Perforated Lower
#0921665	#2 Regular Perforated Lower
#0921477	#3 Regular Perforated Lower
#0921478	#4 Regular Perforated Lower
#0921666	#5 Regular Perforated Lower
#0921667	Set 106 (10) 5 Upper, 5 Lower

#0921186	#0 Regular Solid Upper (XL)
#0921239	#0 Regular Solid Lower (XL)

#0921479	#0 Regular Perforated Upper (XL)
#0921180	#0 Regular Perforated Lower (XL)

#0921720	#15 Pedo/Ortho Solid Upper
#0921721	#16 Pedo/Ortho Solid Upper
#0921722	#17 Pedo/Ortho Solid Upper
#0921723	#18 Pedo/Ortho Solid Upper
#0921724	#15 Pedo/Ortho Solid Lower
#0921725	#16 Pedo/Ortho Solid Lower
#0921726	#17 Pedo/Ortho Solid Lower
#0921727	#18 Pedo/Ortho Solid Lower
#0921719	Set 107 (8) 4 Upper, 4 Lower

#0921872	#15 Pedo/Ortho Perforated Upper
#0921873	#16 Pedo/Ortho Perforated Upper
#0921874	#17 Pedo/Ortho Perforated Upper
#0921875	#18 Pedo/Ortho Perforated Upper
#0921876	#15 Pedo/Ortho Perforated Lower
#0921877	#16 Pedo/Ortho Perforated Lower
#0921878	#17 Pedo/Ortho Perforated Lower
#0921879	#18 Pedo/Ortho Perforated Lower
#0921718	Set 108 (8) 4 Upper, 4 Lower

#0921447	#30 Partial Solid Upper Left and Lower Right
#0921448	#31 Partial Solid Upper Right and Lower Left
#0921449	#32 Partial Solid Anterior Lower
#0921450	#33 Partial Solid Anterior Upper
#0921451	Set 110 (4) #s 30, 31, 32, 33

#0921452	#30 Partial Perforated Upper Left and Lower Right
#0921453	#31 Partial Perforated Upper Right and Lower Left
#0921454	#32 Partial Perf. Anterior Lower
#0921455	#33 Partial Perf. Anterior Upper
#0921456	Set 109 (4) #s 30, 31, 32, 33

#0921866	Twist Tray Perforated
#0921868	Twist Tray Solid

IMPRESSION MATERIALS



PLASTOPASTE™

Plastopaste is a zinc oxide eugenol impression paste (light body) used for preliminary compound impressions, mucostatic impressions, rebasing dentures, and denture relines. Plastopaste easily mixes in 45 seconds to a consistent, extra creamy mixture, and sets hard in 2 minutes 30 seconds without distortion. Plastopaste is non-sticky, wipes easily from patient's face and hands, and can be repositioned without any visible cracks or lines. This formula has a pleasant taste and causes no eugenol burn.

Mixing Ratio: 13 g base to 5.6 g catalyst. I To disinfect: Immerse in a 1:213 iodophor for 10 minutes or 2% glutaraldehyde solution for 10 hours, prepared according to the manufacturer's instructions.

Tips: Plastopaste can be used as a temporary cement.

Pâte à empreinte à l'eugénol et à l'oxyde de zinc (fluide) • Zinkoxid-Eugenol-Abdruckpaste (flüssig) • Pasta per impronte in ossido di zinco eugenolo (molto fluida) • Pasta de impressão de óxido de zinco/eugenol (consistência leve) • Pasta de impresión de eugenol de óxido de zinc (cuerpo ligero)



COMPONENTS		
#0921590	Standard Kit	case: 47
<div><div>7.5oz. (212g) Base</div><div>3.25oz. (92g) Catalyst</div></div>		



H Ships Hazardous

RIMSEAL™

Rimseal is a PEMA acrylic peripheral impression material and bite registration material that can be used with custom or plastic impression trays for full denture impressions and border moulding (muscle trim). Rimseal extends to the proper peripheral height without slumping. This soft material records peripheral detail without distortion and sets to become a rigid part of the tray. Rimseal is yellow in colour, providing a distinct contrast colour from impression materials.

COMPONENTS		
#0921710	Standard Kit (Yellow Powder)	
<div><div>8 oz. (227g) Powder</div><div>6 oz. (117ml) Liquid</div></div>		





H Ships Hazardous

SAPPHIRE™

Sapphire is a PEMA acrylic peripheral impression material and bite registration material that can be used with custom or plastic impression trays for full denture impressions and border moulding (muscle trim). Sapphire extends to the proper height without slumping. This soft material records peripheral detail without distortion and sets to become a rigid part of the tray. No-burn, no-sting formula ensures patient acceptance.

Bulk packaging available.

COMPONENTS		
#0921740	Standard Kit (Blue Powder)	
<div><div>8 oz. (227g) Powder</div><div>6 oz. (117ml) Liquid</div></div>		



SUPERGEL®, SUPERGEL® FRESH & SUPERGEL® MAGIC

Supergel is a premium blended alginate material that mixes readily to a smooth, creamy consistency. Its slow dehydration and low surface tension provide consistent, accurate, and bubble-free detailed reproductions. Pleasing sassafras flavour, controlled flow, heavy body, and gag reflex reduction provide greater patient acceptance. **Supergel Fresh** has all the benefits of our traditional alginate but is a dust-controlled formula for less airborne particles when opening and mixing. Supergel Fresh is blue in colour and berry-cherry flavoured. **Supergel Magic** is a premium blended alginate that mixes readily to a smooth, creamy consistency. Purple during mixing changing to white when set. Sassafras flavour.

Mixing Ratios: 9 g powder + 25 ml water = partial impression. 18 g powder + 50 ml water = full impression.

Poudre d'Empreinte a Base d'Alginate • Polvo Para Molde de Alginato • Alginat-Abdruckmaterial

SUPERGEL		SUPERGEL FRESH (DUST-CONTROLLED FORMULA)		SUPERGEL MAGIC (CHROMATIC FORMULA)	
#0921825	1 lb (454 g) Can - Regular Set	#0921822	1 lb (454 g) Can - Regular Set	#0921812	1lb Can - Regular Set
#0921826	1 lb (454 g) Can - Fast Set	#0921823	1 lb (454 g) Can - Fast Set	#0921813	1lb Can - Fast Set
#0921818	1 lb (454 g) Pouch - Regular Set	#0921827	1 lb (454 g) Pouch - Regular Set		
#0921819	1 lb (454 g) Pouch - Fast Set	#0921830	1 lb (454 g) Pouch - Fast Set		

SUPERPASTE™

Superpaste is a heavy body zinc oxide eugenol impression paste used for taking soft tissue impressions and final impressions with a denture. Superpaste combines the essentials of controlled flowability and “no sag” body to give unexcelled reproduction of the basal tissues. Its non-sticky formula sets hard in approximately 2 minutes to prevent bending or distortion and has no taste or burning sensation for patient comfort.

Tips: Superpaste is an good alternative for taking fully extended peripheral impressions on persons with a sensitivity to acrylics.







Pâte à empreinte à l'eugénol et à l'oxyde de zinc (corps zourd) • Zinkoxid-Eugenol-Abdruckpaste (schwerer Körper) • Pasta per impronte in ossido di zinco eugenolo (il corpo pesante) • Pasta de impressão de óxido de zinco/eugenol (corpo pesado) • Pasta de impresión de eugenol de óxido de zinc (del cuerpo pesado)

COMPONENTS		
#0921850	Standard Kit	case: 47
<div><div>7.5oz. (212g) Base</div><div>3.25oz. (92g) Catalyst</div></div>		
#0921851	7.5 oz (212 g) Tube Base	
#0921852	3.25 oz (92 g) Tube Catalyst	



SUPERSIL™

Supersil is a polyvinylsiloxane impression material recommended for crown and bridge impressions, edentulous impressions, or other impressions requiring fine detail and high stability. Supersil is an addition silicone type material that is odourless and tasteless. Using a 1:1 automix ratio, Supersil dramatically reduces working time and allows for easy clean-up. Supersil allows for multiple pours with no shrinkage and will remain dimensionally stable for several weeks. Available in light body (yellow) or heavy body (pink).

COMPONENTS		
#0921040	Light Body Tube Kit (Yellow)	
#0921043	Heavy Body Tube Kit (Pink)	
<div><div>95g Base</div><div>95g Accelerator</div><div>.3 oz (9ml) Adhesive</div><div>Mixing Pad</div></div>		
#0921046	.3 oz (9 ml) Bottle Adhesive	
COMPONENTS		
#0921055	Light Body Cartridge Kit (Yellow)	
#0921037	Heavy Body Cartridge Kit (Pink)	
<div><div>2 - 50ml Cartridges</div><div>Mixing Tips</div></div>		
#0921414	Static Mixing Tips (48) 5 inch cartridges	



PLASTER

PLASTOGUM®

Plastogum is high-grade modified durable dental plaster. Plastogum can be used for full corrective denture impressions, wash impressions using a custom tray, preliminary compound impressions, or impressions for rebasing dentures. When mixed with water, it forms a putty-like consistency paste. Plastogum does not displace soft tissue, shrink after setting, flake, or soften and crumble from moisture absorption. Controlled setting time enables ease in border moulding. Pleasant flavour.

#0921580	2 lb (908 g)
----------	--------------



ORTHODONTIC

ACRYLIC RESINS

 **Neocryl™**
Orthodontic Acrylic Resin

 **Neocryl™ Splint**
Orthodontic Acrylic Resin
for Splint Fabrication

ORTHODONTIC ACCESSORIES

 **Bo-Boxes™**
Orthodontic Retainer Cases

 **OrthoGlo™**
Glitter

POLISHES

 **Glaze™**
Self-cure Resin Polish

 **Gorilla Seal™**
Ultra-clear Light-cure Surface Sealant

ACRYLIC RESINS




 Ships Hazardous

NEOCRYL™

Neocryl is designed to create brightly coloured appliances that encourage patient compliance. If coloured liquid is used, the appliance will be translucent. If coloured powder is used, the appliance will be opaque. Neocryl powders are formulated extra fine for quick absorption, creating a stronger, more durable appliance, and eliminating slumping or running of the material during fabrication. Longer working time gives the technician greater design and detail time. All Neocryl products are non-yellowing.

Recommended for use with: OrthoGlo™ | Bulk packaging available.
Tips: Neocryl liquids and powders can be blended to achieve endless possibilities.

#0921610 Neocryl Powder Kit




- | | |
|-----------------|----------------|
| Go Green* | Poppin' Pink* |
| Real Red* | Winter White |
| Boldly Blue* | Night Black |
| Yelling Yellow* | Crystal Clear* |

Neocryl Powder Shade Guide
*INCLUDED IN POWDER KIT

BULK POWDER

#0921612	1 lb (454 g) Powder, Crystal Clear
#0921613	1 lb (454 g) Powder, Go Green
#0921614	1 lb (454 g) Powder, Real Red
#0921615	1 lb (454 g) Powder, Boldly Blue
#0921616	1 lb (454 g) Powder, Yelling Yellow
#0921617	1 lb (454 g) Powder, Poppin' Pink
#0921645	1 lb (454 g) Powder, Winter White
#0921646	1 lb (454 g) Powder, Night Black

#0921611 Neocryl Liquid Kit



- | | |
|--------------------|---------------|
| Candy Red* | Very Violet |
| Berry Blue* | Denture Pink |
| Sunny Yellow* | Crystal Clear |
| Grasshopper Green* | |

Neocryl Liquid Shade Guide
*INCLUDED IN LIQUID KIT

BULK LIQUID

#0921624	8 oz (236 ml) Liquid, Crystal Clear
#0921620	8 oz (236 ml) Liquid, Grasshopper Green
#0921621	8 oz (236 ml) Liquid, Candy Red
#0921622	8 oz (236 ml) Liquid, Berry Blue
#0921623	8 oz (236 ml) Liquid, Sunny Yellow
#0921644	8 oz (236 ml) Liquid, Very Violet
#0921640	8 oz (236 ml) Liquid, Denture Pink




 Ships Hazardous

NEOCRYL™ SPLINT

Extra strong, high density acrylic resin designed to make orthodontic splints that help patients who have TMJ or a tendency to grind at night. It will not chip, break, or wear easily and is crystal clear in the patient's mouth. Better liquid absorption prevents running and slumping, resulting in a stronger, more durable appliance that will not yellow.

Bulk packaging available.

#0921636 Neocryl Splint Standard Kit



ORTHODONTIC ACCESSORIES

BO-BOXES™

Bo-Boxes are sturdy, rigid orthodontic retainer cases for storing and protecting appliances. Bo-Boxes are available in five bright colours and provide a pre-moulded identification area for the patient's name and phone number. Ventilation holes allow air to circulate to keep the appliances dry. All cases are medium-sized.

Outside dimensions: 3.125 in wide x 2.75 in tall x 1.375 in deep.

#0921518	Blue Cases (10)
#0921519	Green Cases (10)
#0921520	Red Cases (10)

#0921521	Orange Cases (10)
#0921522	Pink Cases (10)

#0921517	Assorted Cases (10)
#0921523	Assorted Cases (100)



ORTHOGLO™

Combine OrthoGlo Glitter with Neocryl Orthodontic Resin materials to make the brightest appliances. Patients enjoy choosing from exciting gold or silver glitter to customize their appliances.

#0921618	1 oz Bottle Gold
#0921619	1 oz Bottle Silver



POLISHES

GLAZE™

Glaze is a fast-setting, self-cure resin which reduces polishing time and dries to a clear, hard surface with a high gloss finish. Glaze can be brushed onto a clean composite or acrylic material to achieve a polished surface that enhances the appliance and reduces the need for additional polishing.

Recommended for use on provisional, relin, and orthodontic appliances. Glaze fills microscopic voids and reduces bacterial growth. Sets in 20 seconds.

Finierlack Für die Orthodontie, Zahnprothesen, Kronen, Unterfütterungen und Provisionen • Glaçure de finissage Pour appareils d'orthodontie, prothèses, couronnes, rebasages et restaurations temporaires
• Smalto per finitura Per apparecchi ortodontici, dentiere, capsule, ribasature e provvisori • Verniz de acabamento Para aparelhos ortodônticos, dentaduras, coroas, revestimentos e provisórias •
Esmalte para acabado Para aparatos ortodônticos, prótesis dentales, coronas, revestimientos y aplicaciones temporales

#0921236 Standard Kit



GORILLA SEAL™


Gorilla Seal is a crystal clear light-cure surface sealant with a smaller oxygen-inhibited layer that allows the sealant to cure to a high gloss or glossy finish for reduced polishing. Aesthetically pleasing in its 100% crystal clear form, this resin agent can seal over veneers, temporary crowns, composite restorations and I.D. or decorative designs. Gorilla Seal may reduce staining and sensitivity in cases such as tooth bleaching, cervical erosion and veneers.

Tips: Use a halogen/LED with halogen capability for curing.

#0921050 Standard Kit (Bottle)



0921051 Standard Kit (Syringe)



4 - 1.2ml syringes
Gorilla Seal



PREVENTIVE

FLUORIDES

- ☐ **Gelato APF Fluoride Gel & Neutral PH Fluoride Gel**
Fluoride Treatment
- ☐ **Gelato APF Foam Fluoride**
Fluoride Treatment
- ☐ **Gelato Home-care Stannous Fluoride Gel & Vivid Neutral Sodium Fluoride Gel**
Home Fluoride Treatment
- ☐ **Gelato Perio Maintenance Rinse**
Fluoride Maintenance

PROPHYLAXIS

- ☐ **Fortis Disposable Propy Angle**
Propy Angle
- ☐ **Gelato Propy Paste**
Propy Paste
- ☐ **ProphyBrite™**
Air Polishing Unit
- ☐ **Prophy Powder**
Sodium Bicarbonate Powder
- ☐ **Prophylaxis Polishing Brushes & Cups**
- ☐ **Nite Bite**
- ☐ **Pit & Fissure Sealant with Fluoride**
Pit & Fissure Sealant
- ☐ **Tongue Scrapers**

FLUORIDE TRAYS

- ☐ **Hinged EZ Fluoride Trays**
Fluoride Trays
- ☐ **Mizzy Prevent Trays**
Fluoride Trays

SEALANTS

- ☐ **Aegis® Pit & Fissure Sealant with ACP**
Pit & Fissure Sealant

GELATO HOME-CARE STANNOUS FLUORIDE GEL & GELATO HOME-CARE VIVID NEUTRAL SODIUM FLUORIDE GEL

Gelato Home-Care Fluorides are a great way for patients to protect their teeth affordably and effectively at home.

The Gelato Home Care Fluoride Brush-on Gel reduces sensitivity and plaque biofilm, while inhibiting the cariogenic micro process. Gelato Home Care can be easily applied with a tooth brush or a fluoride tray. Contains Xylitol, which is a proven inhibitor of plaque biofilm formation. Available in 0.40% Stannous Fluoride Brush-on Gel or 1.1% Neutral Sodium Fluoride Brush-on Gel.

0.40% STANNOUS FLUORIDE BRUSH-ON		1.1% NEUTRAL SODIUM FLUORIDE BRUSH-ON	
#24-10477	4.3 oz. (.12 L) Mint	#24-09177	4.3 oz. (.12 L) Cherry
#24-10577	4.3 oz. (.12 L) Red Berry	#24-09277	4.3 oz. (.12 L) Mint

FLUORIDE TRAYS

HINGED EZ FLUORIDE TRAYS

These Hinged EZ Fluoride Trays are an excellent way to dispense fluoride treatments with upper and lower arches attached together. The tray is made of closed celled, 4 lb. density foam for good rigidity. With excellent anatomical depth for gingival line, the tray promotes full coverage of the fluoride applications.

#9525720	Small - Yellow (100)	#9525390	Medium - White (100)
#9525890	Small/Medium - Yellow/White (100)	#9525050	Large - Blue (100)

MIZZY PREVENT TRAYS

Mizzy Prevent Trays include exclusive 'fold-away' hinge that prevents gag reflex and assists in easy separation for single-tray use. The tray also has extra-long locking tabs to ensure cleanliness and ease of insertion and removal. The anatomical design and buccal ridge promotes fluoride coverage of the entire tooth while contoured walls reduce overflow and prevent ingestion of fluoride.

Exclusive 'combo' assortment provides all three sizes in one convenient, re-sealable, freestanding bag. The generous 100-count bag offers a substantial saving over the leading national brands.

#7000260	Dual, Small (100)	#7000262	Dual, Large (100)
#7000261	Dual, Medium (100)		

PROPHYLAXIS

FORTIS DISPOSABLE PROPHY ANGLE

The Fortis Propy Angle has been engineered to run longer and smoother with less chatter, providing superior comfort with less hand fatigue than traditional angles. Independent testing has confirmed that the design of the Fortis angle out-performed all other disposable prophy angles with respect to ease of use.

These prophy angles have a second generation gear design with ergonomic grips that reduce the potential fatigue through the hand and arm. The cups are 100% latex free and are available in soft cup, firm cup, pedo cup and tapered brush.

#24-13077	Firm Cup (100)	#24-13277	Pedo Cup (100)
#24-13177	Soft Cup (100)	#24-13377	Tapered Brush (100)

PROPHYLAXIS POLISHING BRUSHES & CUPS

These Propy Brushes & Cups are used for polishing and removal of heavy stains. The brushes are made of fine nylon bristles assembled to prevent the possibility of bristles pulling out. The cups are 100% latex free and made of FDA-listed synthetic material. The cups are webbed to conform to tooth surface and retain polishing material. Available in snap-on, latch type and screw-type.

#24-40014	Propy Brush RA, Flat (144)	#24-40614	Propy Cups Screw Type (144)
#24-40514	Propy Brush RA, Latch Type (144)	#24-40714	Propy Cups Snap-on (144)

FLUORIDES



CE 0120

GELATO APF FLUORIDE GEL & NEUTRAL PH FLUORIDE GEL

Gelato APF Gels are an economical and quick application fluoride treatment. The gel is a 60-second acidulated phosphate fluoride that contains 1.23% fluoride ion. The low PH provides maximum fluoride uptake in less than 60 seconds and the improved thixotropic properties ensure closer adherence to enamel and interproximal areas. The smooth creamy gel will not run and thickens during treatment to prevent patient gagging. The gel is gluten free and offered in many different flavours for patient satisfaction and there's no bitter aftertaste following application. Flip Top Bottle.

Neutral PH formula fluoride gels provide total treatment for a wide patient profile. It is ideal for patients with cosmetic restorations, porcelain or resin. It provides the perfect non-acidic balance, allowing patients to tolerate and accept fluoride applications with comfort with a 4 minute thixotropic formula.

		NEUTRAL PH FLUORIDE GELS	
#24-08677	16 oz. (.48 L) Bubble Gum	#24-01777	16 oz. (.48 L) Bubble Gum
#24-03377	16 oz. (.48 L) Strawberry	#24-01677	16 oz. (.48 L) Strawberry
#24-08277	16 oz. (.48 L) Cherry	#24-00777	16 oz. (.48 L) Dye-Free Mint
#24-10177	16 oz. (.48 L) Cotton Candy		
#24-01877	16 oz. (.48 L) Grape		
#24-10677	16 oz. (.48 L) Mango Smoothie		



GELATO APF FOAM FLUORIDE

Gelato Foam Fluoride is a 1.23% Acidulated Phosphate Fluoride formula that provides excellent, consistent and efficacious coverage. The light foam consistence makes for easy, fast and more pleasant application. The fluoride stays in under bite pressure to eliminate patient gagging.

#24-08977	4.4 oz. (.13 L) Bubble Gum	#24-05577	4.4 oz. (.13 L) Strawberry
#24-08877	4.4 oz. (.13 L) Cherry	#24-04177	4.4 oz. (.13 L) Grape

GELATO PERIO MAINTENANCE RINSE

The 0.63% Stannous Fluoride Perio Rinse has been shown to be an effective therapy to prevent tooth decay while strengthening and protecting roots. The antimicrobial and alcohol-free Gelato Perio Rinse relieves tooth sensitivity, reduces gingival inflammation, helps inhibit plaque build-up, prevents demineralization and promotes re-mineralization.

#24-11877	10 oz. (.29 L) Mint	#24-10377	10 oz. (.29 L) Raspberry
-----------	---------------------	-----------	--------------------------





Disposable
Cups



Jars



WATCH THE
TECHNIQUE VIDEO!



GELATO PROPHYLAXIS PASTE

Gelato has been a forerunner of great-tasting prophylaxis paste, bringing on new flavours such as Pina Colada and Orange Sherbet but what also makes Gelato Prophylaxis Cups great is a non-splatter formula that contains 1.23% APF with a unique blend of cleaning and polishing agents.

Gelato paste comes in boxes of 200 individual disposable cups for convenient use. The disposable cups also eliminate cross contamination and include a proph ring for ease of application on the patient's teeth. For more options, the paste also comes in 6 ounce and 12 ounce jars.

Gelato Prophylaxis Paste is offered in 4 different grits: fine, medium, coarse and x-coarse. The fine grit is best for light stain removal and amalgam polishing. It's recommended to use the fine grit on children. For normal cleaning and polishing the medium girt is an ideal choice. The coarse grit is used in medium to heavy stain and plaque removal.

JARS

#24-10877	12 oz (.35 L)	Fine Bubble Gum
#24-09077	12 oz (.35 L)	Medium Cherry
#24-08177	12 oz (.35 L)	Medium Mint

#24-08077	12 oz (.35 L)	Medium Piña Colada
#24-03977	12 oz (.35 L)	Coarse Cherry
#24-04077	12 oz (.35 L)	Coarse Mint

CUPS

#24-02877	Fine Bubble Gum (200)
#24-02377	Medium Bubble Gum (200)
#24-02477	Medium Cherry (200)
#24-02277	Medium Mint (200)
#24-02577	Medium Piña Colada (200)
#24-11677	Medium Orange Sherbert (200)
#24-15877	Medium Raspberry (200)
#24-02677	Medium Assorted* (200)
#24-07577	Coarse Bubble Gum (200)

#24-02977	Coarse Cherry (200)
#24-03077	Coarse Mint (200)
#24-03177	Coarse Piña Colada (200)
#24-11777	Coarse Orange Sherbert (200)
#24-15977	Coarse Raspberry (200)
#24-02777	Coarse Assorted* (200)
#24-05077	X-Coarse Cherry (200)
#24-05177	X-Coarse Mint (200)
#24-15177	Prophy Ring (200)

*Assorted flavours: Bubble Gum, Cherry, Mint and Piña Colada

PROPHYBRITE™

The ProphylBrite air polisher provides quick and gentle patient care. This unit delivers a fine sodium bicarbonate and water slurry for easy removal of plaque and stains. Ideal for cleaning and preparing teeth prior to placement of orthodontic brackets, sealants, bleaching, and fluoride treatments. Offers chair side convenience. Unit is easy to install, light weight, and portable from operator to operator. Connects to most four-hole handpiece coupler/air and water connection. A specific four-hole hose connection should be designated for optimal results. Comes with a DVD.

Tips: Nozzles are autoclavable. Do not autoclave nozzles if clogged. Adaptors are for use with Bosworth ProphylBrite only.

“The Bosworth ProphylBrite™ is both ergonomic and efficient. The light weight hand piece is portable and it can be transferred from one operator to another with no special water source hookups. The Bosworth ProphylBrite is not only easy to fit into your practice, but easy to use! Dr. Sheri B. Doniger”

#16670R	Air Polishing Unit Kit (Red)
#16670B	Air Polishing Unit Kit (Black)
#16670LB	Air Polishing Unit Kit (Blue)
#16670N	Air Polishing Unit Kit (Orange)

(includes Nozzle (2), Air Polisher, Powder Chamber Cap (2), Cleaning Files (2.5 cm, 8 cm, and 10 cm), O-Rings (cap, powder and water), Bosworth Prophy Powder (3 oz, 85 g sample)

Dimensions: 7.5 in x 2 in x 2.4 in (19 cm x 5 cm x 6 cm)
Water Delivery: 40-60ml/min • Service Pressure: 65-70 p.s.i.

ADAPTORS	
#16678	Adaptor - Kavo
#16684	Adaptor - W&H-Adec
#16680	Adaptor - NSK Mach-QD
#16683	Adaptor - NSK Ti-Max
#16671	Adaptor - Bien Air
#16686	Adaptor - DentalEZ

COMPONENTS

#16677	Nozzle
#16681K	Spare Parts Set (1 Powder Chamber Cap & 1 set of O-rings)

PROPHY POWDER

Less abrasive than prophyl paste, use sodium bicarbonate prophyl powder on enamel to remove stubborn stains, plaque and soft debris. Bosworth's prophyl powder is compatible with most air polishing units.

#16690	10 oz Grape	#16692	10 oz Orange	#16631	10 oz Raspberry
#16691	10 oz Spearmint	#16630	10 oz Strawberry	#16632	10 oz Lemon-Lime

#16696	Tri-Pack (3 oz Grape, 3 oz Orange, 3 oz Spearmint)
#16636	Tri-Pack (3 oz Lemon-Lime, 3 oz Raspberry, 3 oz Strawberry)

SEALANTS


AEGIS® PIT & FISSURE SEALANT WITH ACP

Aegis Pit and Fissure Sealant contains the “smart material” Amorphous Calcium Phosphate (ACP). The ACP filler is resilient and flexible, creating a stronger and longer lasting sealant. Aegis P&F has low shrinkage, good margin, and a slow rate of stress development so as to not chip off. Aegis P&F can be used with any light-curing device. Aegis' controlled flowability keeps the sealant on the tooth while completely filling the pits and fissures. When properly maintained and placed, Aegis forms a chemical barrier that may help keep patients free from carious lesions on the occlusal surfaces. Sealant cures to opaque on enamel.


Recommended for use with: Acid Etch and Bosworth ProphylBrite®.

Vernis dentaire photopolymérisable avec phosphate de calcium amorphe (ACP) • Lichtaushärtende Versiegelung mit ACP • Sigillante fotopolimerizzabile con ACP (fosfato di calcio amorfo) • Resina selante de cura por luz con ACP • Sellador de curado por luz con ACP (fosfato de calcio amorfo)


#0921348 Standard Kit



#0921327 Clinic Kit



#0921328 Refill Kit




PIT & FISSURE SEALANT WITH FLUORIDE

Pit & Fissure Sealant is a light-cured resin with Fluoride. The controllable flowability keeps the sealant on the tooth while completely filling the pits and fissures. This sealant cures to an opaque film on enamel.

Recommended for use with: Acid Etch.

#092156 Standard Kit



TMD APPLIANCE

NITE BITE

For the first time ever you can now stop bruxism, clenching, and jaw joint discomfort (TMD). NiteBite's® patented technology utilizes the proprioceptive senses to initiate the jaw opening reflex within the Freeway Space, allowing the jaw to relax at the physiologic rest position. Thus, bruxism stops and its related symptoms are relieved.

NiteBite is designed to fit almost every patient. It can be fitted in only five minutes. No impression is required and no second patient visit is needed. So treating your patients with the Nitebite is quick, easy and cost effective. In fact, Nitebite can and should be used as the first appliance you choose to treat any patient who has symptoms of TMD.

#9541130 NiteBite Kit (5)

TONGUE SCRAPERS

TONGUE SCRAPERS

Did you know that more than 85% of bad breath originates in the mouth, not the stomach? Tongue Scrapers provide a handy way to clean the tongue's surface and make sure that bad breath bacteria are scraped away. Cleaning on a regular basis is known to significantly reduce bacteria and food particles on the tongue. Reduces mouth odours, overall bacteria and helps to keep the mouth fresh and clean.

#9541130	Assorted (100)
#9541140	Light Blue (100)

#9541160	Dark Blue (100)
#9541190	Green (100)

#9541220	Red (100)
#9541240	White (100)



WATCH THE
TECHNIQUE VIDEO!



PROVISIONAL

TEMPORARY RESIN ACRYLICS

- ☐ **Trim® & Trim® II**
PEMA Temporary Resin Acrylic
- ☐ **Trim® Plus**
PMMA Temporary Resin Acrylic

- ☐ **Ultra Trim®**
Bis-Acryl Temporary Crown
and Bridge Material

POLISHES

- ☐ **Glaze™**
Self-cure Resin Polish
- ☐ **Gorilla Seal™**
Ultra-clear Light-cure Surface Sealant

TEMPORARY RESIN ACRYLICS



WATCH THE
TECHNIQUE VIDEO!



PART 1



PART 2



WATCH THE
TECHNIQUE VIDEO!



PART 1



PART 2



TRIM® AND TRIM® II

Trim and Trim II are fast-setting, non-irritating PEMA temporary resin acrylic crown and bridge materials formulated for extended use. The powder and liquid in each kit are identical in formula and the powders can be mixed to obtain speciality colours. The resin is dimensionally stable, will not distort when withdrawn from undercuts while setting and is suitable for use as a bite registration material. Trim is available in two shades (light and dark). Trim II is available in six shades for greater diversity in achieving the perfect shade match for patients.


Recommended for use with: Glaze™, Gorilla Seal™, Masque™ and Pressure Pot I Bulk packaging available.
Mixing Ratio: 6.3 g of powder to 3 ml of liquid. This amount would make 9 full jacket crowns. I Working time: 3-4 minutes.

Tips: Order Trim VW (Very White) powder in addition to the standard kit to produce a temporary as white as bleached teeth.


Note: Trim and Trim II cannot be mixed with Trim Plus materials.

Kaltpolymerisat für Provisorische Kronen und Brücken • Resina para Caronas y Puentes Temporales • Resina pour Couronnes et Bridges Provisoires • Resine per Ponti e Corone Provvisorie

#0921900 Trim Standard Kit case: 20




#0921090 Trim II Standard Kit case: 15



#0921905 Trim 1 lb (454 g) Bulk Kit, All Light Powder

#0921906 Trim 1 lb (454 g) Bulk Kit, All Dark Powder



COMPONENTS

- #0921901 Trim Light Powder 3 oz (85 g)
- #0921902 Trim Dark Powder 3 oz (85 g)
- #0921999 Trim VW (Very White) Powder 3 oz (85 g)
- #0921092 Trim II 1.5 oz (42 g) Powder, Shade B1 (59)
- #0921105-60 Trim II 1.5 oz (42 g) Powder, Shade A1 (60)
- #0921105-66 Trim II 1.5 oz (42 g) Powder, Shade A3 (66)
- #0921093 Trim II 1.5 oz (42 g) Powder, Shade B2 (61)
- #0921094 Trim II 1.5 oz (42 g) Powder, Shade A2 (62)
- #0921095 Trim II 1.5 oz (42 g) Powder, Shade D3 (65)
- #0921096 Trim II 1.5 oz (42 g) Powder, Shade C2 (69)
- #0921097 Trim II 1.5 oz (42 g) Powder, Shade B4 (77)
- #0921907 Trim 1 lb (454 g) Light Powder
- #0921908 Trim 1 lb (454 g) Dark Powder
- #0921100 Trim II 1 lb (454 g) Bulk Powder Select Shade
- #0921091 4 oz (118 ml) Liquid
- #0921903 6 oz (177 ml) Liquid
- #0921903s 6 oz (177 ml) Fast Liquid
- #0921909 8 oz (236 ml) Liquid
- #0921914 32 oz (946 ml) Liquid
- #0921089 3 ml Pipette Dropper (15)

TRIM® PLUS


Trim Plus is a PMMA chairside temporary resin acrylic for use in fabricating temporary crowns, bridges, partials, and splints. This PMMA formula is colour-stable for extended life and is available in shades B1 (59), B2 (61), A1 (60), A2 (62), D3 (65), A3 (66), B3 (67), 68, C2 (69), B4 (77), A3.5 (81), and A4 (87) for greater selection. Shades provide a natural appearance when matching patient requirements.

Recommended for use with: Glaze™, Gorilla Seal™, Masque™ and Pressure Pot I Bulk packaging available.

Note: Trim Plus cannot be mixed with Trim and Trim II materials.

Kaltpolymerisat für Provisorische Kronen und Brücken • Resina para Caronas y Puentes Temporales • Resina pour Couronnes et Bridges Provisoires • Resine per Ponti e Corone Provvisorie

#0921930 Standard Kit



COMPONENTS

- #0921931 1.5 oz (42 g) Powder Select Shade
- #0921931s 1.5 oz (42 g) Fast Set Powder Select Shade
- #0921936 4 oz (118 ml) Liquid

BULK


- #0921933 1 lb (454 g) Bulk Powder Select Shade
- #0921933s 1 lb (454 g) Bulk Fast Set Powder Select Shade
- #0921937 8 oz (236 ml) Liquid

ULTRA TRIM™

Ultra Trim is an auto-mixed, cartridge dispensed bis-acryl temporary crown and bridge material that will greatly reduce in-office procedure time. Ultra Trim's 10:1 base to catalyst mixing ratio provides excellent strength and greater flexibility without brittleness. (The higher the base content, the stronger the restoration). Ultra Trim will bond to itself if additional material or repair is necessary. Ultra Trim's non-irritating formula is free of any harsh monomers or offensive odours. It can be used in virtually any matrix with minimal shrinkage for excellent fit and marginal adaption. Ultra Trim sets in the mouth in approximately 2 minutes and completely cures in 5 minutes. Ultra Trim is colour-stable with perfect translucency.

Recommended for use with: E-Z Squeeze™ Universal Dispenser, Glaze™, Gorilla Seal™, Masque™ and Pressure Pot

Note: Cartridge should be bled before every use to ensure base and catalyst go into the dispensing tip evenly. A 4:1 tip and dispenser must be used with the 31 g cartridges.



COMPONENTS

- #0921803 Starter Kit A1 (Light)
- #0921804 Starter Kit A2 (Universal)
- #0921805 Starter Kit A3.5 (Dark)
- #0921806 76 g A1 Cartridge (Light)
- #0921807 76 g A2 Cartridge (Universal)
- #0921808 76 g A3.5 Cartridge (Dark)
- #0921699 15 - 10:1 Mixing Tips
- #0921236 Glaze Standard Kit
- #0921809 10:1 E-Z Squeeze Universal Dispenser

New Hue & VITA Shade Reference Chart

	Trim	Trim II	Trim Plus	Ultra Trim
VW.....B0	•			
59.....B1		•	•	
60.....A1		•	•	•
61.....B2	•	•	•	
62.....A2		•	•	•
65.....D3	•	•	•	
66.....A3		•	•	
67.....B3			•	

	Trim	Trim II	Trim Plus	Ultra Trim
68.....			•	
69.....C2		•	•	
77.....B4		•	•	
81.....A3.5			•	•
87.....A4			•	
*Clear	•			
*Incisal ...			•	

*Made available for characterization and for putting an incisal edge on temporaries.

POLISHES

GLAZE™

Glaze is a fast-setting, self-cure resin which reduces polishing time and dries to a clear, hard surface with a high gloss finish. Glaze can be brushed onto a clean composite or acrylic material to achieve a polished surface that enhances the appliance and reduces the need for additional polishing. Recommended for use on provisional, relene, and orthodontic appliances. Glaze fills microscopic voids and reduces bacterial growth. Sets in 20 seconds.

Finierlack Für die Orthodontie, Zahnprothesen, Kronen, Unterfütterungen und Provisorien • Glaçure de finissage Pour appareils d'orthodontie, prothèses, couronnes, rebasages et restaurations temporaires • Smalto per finitura Per apparecchi ortodontici, dentiere, capsule, ribasature e provvisori • Verniz de acabamento Para aparelhos ortodônticos, dentaduras, coroas, revestimentos e provisórios • Esmalte para acabado Para aparatos ortodônticos, prótesis dentales, coronas, revestimientos y aplicaciones temporales

#0921236 Standard Kit



GORILLA SEAL™


Gorilla Seal is a crystal clear light-cure surface sealant with a smaller oxygen-inhibited layer that allows the sealant to cure to a high gloss or glossy finish for reduced polishing. Aesthetically pleasing in its 100% crystal clear form, this resin agent can seal over veneers, temporary crowns, composite restorations and I.D. or decorative designs. Gorilla Seal may reduce staining and sensitivity in cases such as tooth bleaching, cervical erosion and veneers.

Tips: Use a halogen/LED with halogen capability for curing.

#0921050 Standard Kit (Bottle)



#0921051 Standard Kit (Syringe)



WATCH THE
TECHNIQUE VIDEO!



RESTORATIVE

CEMENTS

 **Aegis® Crown & Bridge**
Crown & Bridge Cement with ACP

 **Masterlink®**
Dual-cure Resin Cement

 **SuperEBA®**
All-purpose Liner & Cement with Alumina

 **Zinc Phosphate Cement**
Permanent Cement

LINERS & VARNISHES

 **Aegis® Liner**
Cavity Liner/Base Material

 **Copaliner®**
Cavity Varnish and Sealant

 **Hydrox®**
Calcium Hydroxide Paste

POLISHES

 **Gorilla Seal™**
Ultra-clear Light-cure Surface Sealant

RESTORATIVE MATERIALS

 **Aegis® V**
Class V Restorative Material

 **Esteem®**
Self-cure Composite Restorative Material

CEMENTS



AEGIS® CROWN & BRIDGE CEMENT WITH ACP



Aegis Crown & Bridge is a syringe dispensed, two paste, dual-cure modified resin polymer material. It is the first crown & bridge resin cement containing an Amorphous Calcium Phosphate (ACP) filler. Aegis C & B is ideal for use as a permanent luting cement for crown & bridge, inlay, and veneer cementation applications, and makes for excellent margin. When exposed to a curing light, Aegis C & B cures in 5-10 seconds.

Ciment à la résine à durcissement double avec de l'ACP • Dualhärtender Kompositzement mit ACP • Cimento in resina a polimerizzazione duale con ACP • Cimento resinoso de dupla polimerização com ACP • Cimento de resina de doble cura con fosfato de calcio amorfo (ACP)

#0921025 Standard Kit (Shade B1 Light)			
			
2 ml syringes Aegis C&B	6ml Aqua Bond™	6ml Dentin Resin	20 Mixing Tip

MASTERLINK®

MasterLink Crown and Bridge is a syringe dispensed, two paste, dual-cure resin cement, ideal for use with gold, ceramic, crowns, inlays and onlays. Masterlink features low viscosity and thin film thickness to provide ideal handling. A non-runny material makes a perfect margin the first time and every time. MasterLink is radiopaque and available in Vita Shade B1 (59), A2 (60), and C1.

#092176	Standard Kit
	
1 - 9 g syringe MasterLink Cement	5 Mixing Tips





SUPEREBA®

SuperEBA is an all-purpose, alumina-fortified material for final cementation of a single crown or inlay, or for use as a base, liner, or a temporary dressing for the relief of pain. SuperEBA's neutral pH and mild formulation provide additional patient comfort during the healing process. Hydrophilic nature and low film thickness provide unsurpassed convenience, making it an excellent choice for patients who are easily irritated. Adheres to tooth structure in the presence of moisture. Not to be used for implants.

Tips: Regular set has a greater than 15 minute slab set, but a shorter in the mouth set of 10-15 minutes. The fast set has a 10-15 minute slab set, but a faster in the mouth set of 4-8 minutes. Film thickness 25 microns.

Bulk packaging available.

ZOE-Unterfütterungsmaterial und -Zement • Revêtement et ciment ZOE universels • Cimento e sottofondo multiuso ZOE • Cimento e revestimento ZOE para todos os fins • Cimento y revestimiento ZOE para todo propósito


#0921007	Standard Kit	Specify (F) for all Fast Powder or (R) for all Regular Powder			COMPONENTS	
					#0921008	15 g Powder Fast Set
					#0921009	15 g Powder Regular Set
					#0921010	6.5 ml Liquid Only

ZINC PHOSPHATE CEMENT

Zinc Phosphate Cement is specially blended to mix smooth and complete. This cement has excellent bonding strength and works well for crowns, inlays, fixed bridgework, orthodontic bands, and temporary fillings. Powder available in yellow or light yellow.

Tips: As the amount of powder incorporated into the mix is increased, the strength increases. Incorporate the powder incrementally. If a large amount of powder is mixed into the liquid all at once, the temperature at the site of the reaction becomes higher, resulting in an increased reaction time and loss of control of the consistency. For orthodontic bands, a heavy consistency will produce greater film thickness, which is important for band retention.


Bulk packaging available.

#092107	Standard Kit (Specify Yellow or Light Yellow)	COMPONENTS	
		#092108	15 ml Liquid
		#092109	29 g Powder (Specify Yellow or Light Yellow)

LINERS & VARNISHES

AEGIS® LINER WITH ACP



Aegis Liner is a light-cured resin with Amorphous Calcium Phosphate (ACP) formulated cavity liner/ base material. Easy to apply and will chemically bond to all adhesive and composites. Available in a convenient dispensing syringe.

#0921730	Standard Kit
	4 - 1.2ml syringes Aegis Cavity Liner with ACP

COPALINER®

Copaliner cavity varnish and sealant protects the dentin and penetrates, fills, and seals the tubuli. It insulates against thermal, galvanic, and electrical shock and prevents hyperemia, death of the pulp, and discolouration. Copaliner provides a strong bond and an excellent marginal seal that will extend the life of the restoration. Its non-chloroform formula is scratch resistant and virtually insoluble in the oral environment. Over time it will develop dentinal tags for improved performance. Copaliner is ether-based with a biologically mild formula, minimizing post-operative sensitivity. CAN ONLY BE USED UNDER AMALGAM.

Vernis à cavités et agent de scellement • Zur Verwendung unter Amalgamfüllungen sowie als Überzug unter Zahnaufbauten • Vernice per cavità dentarie e sigillante • Barniz y sellador para caries





#0921525	Standard Kit	COMPONENTS	
		#0921526	.5 oz (14 ml) Varnish
		#0921527	.5 oz (14 ml) Solvent

HYDROX®

Hydrox radiopaque calcium hydroxide paste is a rigid, self-setting cavity liner and pulp capping agent. Hydrox eliminates all chemical and thermal sensitivity by sealing dentinal tubules and protecting the pulp. Hydrox serves as a protective barrier that insulates the dentin and pulp from acid-containing cements and restorative materials. Works well under glass ionomers and will not affect amalgam or the polymerization of composites. The superior flowability of Hydrox permits uniform coverage and controlled placement into preparations.

Tips: Hydrox can be used as a temporary cement. Hydrox's high pH level acts to inhibit bacterial growth at the cavity wall. This property is particularly important for the treatment of teeth with complicated crown fractures, open apices, or root perforations.

Pâte à L'hydroxyde de calcium • Calciumhydroxide-Paste • Calcium Hydroxide Pasta • Pasta de Hidróxido Cálcico

#0921071	Standard Kit, Dentin	#0921073	6 Pack Kit, Ivory
#0921072	Standard Kit, Ivory	#0921074	6 Pack Kit, Dentin
			
15g Base	13g Catalyst	6 - 15g Base	6 - 13g Catalyst



POLISHES



GORILLA SEAL™


Gorilla Seal is a crystal clear light-cure surface sealant with a smaller oxygen-inhibited layer that allows the sealant to cure to a high gloss or glossy finish for reduced polishing. Aesthetically pleasing in its 100% crystal clear form, this resin agent can seal over veneers, temporary crowns, composite restorations and I.D. or decorative designs. Gorilla Seal may reduce staining and sensitivity in cases such as tooth bleaching, cervical erosion and veneers.

Tips: Use a halogen/LED with halogen capability for curing.

#0921050 Standard Kit (Bottle)



#0921051 Standard Kit (Syringe)



RESTORATIVE MATERIALS




AEGIS® V WITH ACP

Aegis V is an Amorphous Calcium Phosphate (ACP) filled restorative that is a light-cured material specifically designed for Class V restorations. Provides controlled placement for increased patient comfort and effortless application.

Recommended for use with: Gorilla Seal™

#0921735 Standard Kit



4 - 1.2ml syringes
Aegis V with ACP



COMPONENTS

#0921736 1.2 ml Syringe A1 (2)
#0921737 1.2 ml Syringe A2 (2)
#0921738 1.2 ml Syringe A3.5 (2)
#0921739 1.2 ml Syringe B1 (2)

ESTEEM®

Esteem is a vita shade (A-2), 100% radiopaque, self-cure composite restorative material. The two-paste system, which includes a universal paste and catalyst paste, forms a durable restoration. The hardened surface can be polished to a natural appearance, almost identical to original tooth anatomy. Esteem is stain resistant and has physical properties that produces exceptional holding strength. It matches the full range of tooth structure with a single translucent shade that blends into the surrounding dentition. It mixes easily, sets quickly, and finishes beautifully. Can be used for post and core build-up.

Recommended for use with: Acid Etch, Hydrox®, and Wonder Wedges®

Tips: To prevent cross-contamination and hardening within the jars, opposite ends of the mixing stick should be used to dip into the universal and catalyst jars. I Note: This product was formerly known as “Bosworth Natural”.

#0921560 Standard Kit



14g
Universal Paste



14g
Catalyst Paste

ACCESSORIES

#0921564 Double-Ended Spatulas (100)



Disinfection of surgical
and rotating instruments
zeta 1 ultra - zeta 2
sporex - zeta 2 enzyme



Surface disinfection of
medical devices
zeta 3 ultra - zeta 3 soft
- zeta 3 wipes - zeta 3
wipes TOTAL



Delicate surface
disinfection of medical
devices
zeta 3 foam - zeta 3
wipes POP-UP



Deterrence of washable
surfaces
zeta 4 wash



Disinfection of aspiration circuit
zeta 5 unit



Deterrence of hands
zeta 6 hydra - zeta 6 drygel



Disinfection of impressions
zeta 7 solution - zeta 7 spray

Packaging

C810002	ZETA 1 ULTRA	125 ml bottle
C810000	ZETA 1 ULTRA	1 litre bottle
C810001	ZETA 1 ULTRA	5 litres container with dosage cap
C810011	ZETA 2 SPOREX	900 g container with dosing spoon
C810012	ZETA 2 ENZYME	1.2 kg container with dosing spoon
C810030	ZETA 3 ULTRA	125 ml bottle
C810021	ZETA 3 ULTRA	750 ml bottle with nozzle trigger spray
C810022	ZETA 3 ULTRA	3 litres container with dosage cap
C810029	ZETA 3 SOFT	125 ml bottle
C810023	ZETA 3 SOFT	750 ml bottle with nozzle trigger spray
C810024	ZETA 3 SOFT	5 litres (2x2.5 litres) container with dosage cap
C810032	ZETA 3 SOFT CLASSIC	125 ml bottle
C810027	ZETA 3 SOFT CLASSIC	750 ml bottle with foamer trigger spray
C810028	ZETA 3 SOFT CLASSIC	5 litres (2x2.5 litres) container with dosage cap
C810060	ZETA 3 WIPES	1 tub, 120 wipes each
C810061	ZETA 3 WIPES	1 refill bag, 120 wipes each
C810063	ZETA 3 WIPES TOTAL	1 tub, 120 wipes each
C810062	ZETA 3 WIPES TOTAL	1 refill bag, 120 wipes each
C810025	ZETA 3 FOAM	750 ml bottle with foamer trigger spray

C810026	ZETA 3 FOAM	3 litres container with dosage cap
C810064	ZETA 3 WIPES POP-UP	100 wipes soft pack
C810037	ZETA 4 WASH	3 litres container
C800061	ZETA 5 UNIT	5 litres container
C810044	ZETA 6 HYDRA	125 ml bottle
C810042	ZETA 6 HYDRA	1 litre bottle with dosage cap
C810043	ZETA 6 HYDRA	5 litres container with dosage cap
C810047	ZETA 6 DRYGEL	80 ml bottle
C810045	ZETA 6 DRYGEL	500 ml bottle
C810046	ZETA 6 DRYGEL	1000 ml bottle
C810052	ZETA 7 SOLUTION	125 ml bottle
C810048	ZETA 7 SOLUTION	1 litre bottle
C810049	ZETA 7 SOLUTION	5 litres container with dosage cap
C810053	ZETA 7 SPRAY	125 ml bottle
C810050	ZETA 7 SPRAY	750 ml bottle with foam cap

Accessories

C800090		1 litre sterilization bath
C800092		3 litres sterilization bath
C810070	ZETA HYGIENE NOZZLE	Dosage cap for 3 and 5 litres container



THERMOFORMING MATERIALS



PROFORM NIGHTGUARD

This laminated material was specifically made for night time bruxers. The layering gives the thermoplastic a more abrasion resistant characteristic that will withstand grinding better than traditional soft sheet materials.

Available in :	127mm Square	120mm Round	125mm Round
3.00mm sold per blank	NGL	NGL-120	NGL-125



CLEAR HD FOR NIGHTGUARDS AND FLUORIDE TRAYS

These EVA blanks are of the same formula as our bleach tray material which is slightly firmer and more resilient than our mouthguard material.

Available in :	127mm Square	120mm Round	125mm Round
2.00mm Pk of 5 Blanks	CHD2	CHD2-120	CHD2-125
2.00mm Pk of 25 Blanks	CHD2-L	CHD2-L-120	CHD2-L-125
3.00mm Pk of 5 Blanks	CHD3	CHD3-120	CHD3-125
3.00mm Pk of 25 Blanks	CHD3-L	CHD3-L-120	CHD3-L-125



BLEACH TRAY MATERIAL

These EVA blanks have been specially formulated for the bleaching technique and are slightly firmer and more resilient than our mouthguard material.

Available in :	127mm Square	120mm Round	125mm Round
1.00mm Pk of 5 Blanks	O1BT1	O1BT1-120R	O1BT1-125R
1.00mm Pk of 25 Blanks	O1BT1-L	O1BT1-L-120R	O1BT1-L-125R
1.50mm Pk of 5 Blanks	O1BT15	O1BT15-120R	O1BT15-125R
1.50mm Pk of 25 Blanks	O1BT15-L	O1BT15-L-120R	O1BT15-L-125R



FOAM BACKED BLEACH TRAY

This Bleaching Laminate is ideally suited for bleach trays, it is only 0.50mm thick once formed over the model. It has a foam liner which absorbs the bleach, thus decreasing the bleaching time for the patient.

Available in :	127mm Square	120mm Round	125mm Round
Pk of 5 Blanks	O1BTFB	O1BTFB-120R	O1BTFB-125R



TRAY AND BASE PLATE MATERIAL

These blanks are made from high-impact styrene and are ideal for making custom impression trays and base plates.

Available in :	127mm Square Only
Pks of 25 Blanks White 3.00 mm	WT3
Pks of 25 Blanks White 4.00 mm	WT4
Pks of 25 Blanks Blue 3.75 mm	BT4
Pks of 25 Blanks Pink 2.00 mm	PT2
Pks of 25 Blanks Red 1.50 mm	VBP
















Coping Material

These polypropylene blanks are excellent for fabricating temporary crowns, bridges and copings. The material gives good definition and will not adhere to acrylics.

Available in :	127mm Square	120mm Round	125mm Round
0.50mm Pk of 50 Blanks	VCC	VCC-R120	VCC-R125
1.00mm Pk of 50 Blanks	VCC1	VCC1-R120	VCC1-R125

MOUTHGUARD ORDER CODE EXAMPLE

Material Type	Thickness mm	Colour	Pk Size	Blank Size	Complete code
MG	3	-	L	SQ	Eg: Clear Mouthguard - 3mm - Square - Large pk of 25 = MG 3 L SQ
CM	4	NB	-	120	Eg: Coloured Mouthguard - 4mm - Navy Blue - 120mm diameter = CM 4 NB 120
MC3	5	R/W/R	-	125	Eg: Multi-Coloured - 5mm - Red/White/Red - 125mm diameter = MC3 5 R/W/R 125
SD	4	FANGS	-	SQ	Eg: Special Design - 4mm - Fangs - Square = SD 4 FANGS SQ
FG	4	FIRE	-	120	Eg: Funguard - 4mm - 120mm diameter = FG 4 FIRE 120
GG	4	GR	-	SQ	Eg: Glitter Guard - 4mm - Green - Square = GG 4 GR SQ

 Glass Clear 1, 1.5, 2, 3, 4, 5mm	 Red 2, 4, 5mm, Ref - R	 White 2, 4, 5mm, Ref - W	 Royal Blue 2, 4, 5mm, Ref - B	 Yellow 2, 4, 5mm, Ref - Y	 Green 2, 4, 5mm, Ref - G
 Black 2, 4, 5mm, Ref - BK	 Orange 4, 5mm, Ref - O	 Fluorescent Pink 4, 5mm, Ref - FP	 Fluorescent Green 4, 5mm, Ref - FG	 Fluorescent Orange 4, 5mm, Ref - FO	 Fluorescent Yellow 4, 5mm, Ref - FY
 Fluorescent Red 4, 5mm, Ref - FR	 Maroon 4mm, Ref - M	 Purple 4, 5mm, Ref - P	 Silver 4mm, Ref - S	 Gold 4mm, Ref - GLD	 Navy 4mm, Ref - NB
 Light Blue 4mm, Ref - LB	 Sky 4mm, Ref - SB	 Racing Green 4mm, Ref - RG	 Baby Pink 4mm, ref - SP	 Lavender 4mm, Ref - LA	 Turquoise 4mm, Ref - T
 Lilac 4mm, Ref - LIL	 Brown 4mm, Ref - BR	 Apple 4mm, Ref - A	 Sandstorm 4mm, Ref - SAND	 Multi Colour Madness 4mm, Ref - MAD	 Cool Zebra 4mm, Ref - ZEB
 Inferno 4mm, Ref - INFER	 Cool Denim 4mm, Ref - DEN	 Aurora Lights 4mm, Ref - ALITE	 Moonshine 4mm, Ref - MOON	 Storm 4mm, ref - STORM	 Patriot 4mm, Ref - PAT
 Rainbow 4mm, Ref - RAIN	 Shamrock 4mm, Ref - SHAM	 Fire Storm 4mm, Ref - FIRE	 Camouflage 4mm, Ref - CAMO	 Flamingo 4mm, Ref - FLAM	 Purple Glitter 4mm, Ref - GP
 Green Glitter 4mm, Ref - GG	 Blue Glitter 4mm, Ref - GB	 Pink Glitter 4mm, Ref - GP	 Yellow Glitter 4mm, Ref - GY	 Red Glitter 4mm, Ref - GR	Please note colours and designs are indicative only, an exact match cannot be guaranteed.



CLEAR SOFT MOUTHGUARD BLANKS

A soft EVA blank ideal for fabricating Mouthguards, Flouride trays, Bleaching trays & Nightguards. Thicknesses available 1.00, 1.50, 2.00, 3.00, 4.00, 5.00mm. Blanks supplied in either 127mm Square or 120mm & 125mm Round. Supplied in packs of 5 or 25 blanks, except 5mm which is supplied in packs of 5 and 20 blanks


Ref - MG



COLOURED MOUTHGUARD BLANKS

A soft EVA blank ideal for fabricating Mouthguards, in a range of vibrant colours. Available in a range of thicknesses see colour chart. Blanks supplied in either 127mm Square or 120mm & 125mm Round

Ref - CM



MULTI-COLOURED STRIPED MOUTHGUARD BLANKS

A soft EVA striped blank supplied in any colour combination within our range. So choose your team colours etc and stand out from the crowd. Blanks available in 127mm square or 120mm & 125mm round in either 4.00 mm & 5.00 mm thicknesses. The colour combinations illustrated below are for display purposes only, any colour can be substituted for another colour or clear providing that colour is available in the thickness required.



2 Stripe Example

Ref - MC2



3 Stripe Example

Ref - MC3



4 Stripe Example

Ref - MC4



5 Stripe Example

Ref - MC5



6 Stripe Example

Ref - MC6



7 Stripe Example

Ref - MC7

SPECIAL DESIGNS MOUTHGUARD BLANKS

A soft EVA blank with a unique special design. Blanks available in 127mm square or 120mm & 125mm round in either 4.00 mm & 5.00 mm thicknesses. The colours of the special designs maybe substituted for an alternative colour within our range.



Ref - SD St Georges



Ref - SD Union Jack



Ref - SD Fangs



Ref - SD Goofy



Ref - SD Checkered Example




Ref - SD Stars and Stripes

Note the stars for the stars and stripes blank are stickers and will require laminating or sealing using our centri-fuse kit

FUNGUARD BLANKS

A soft 4.00mm EVA blank ideal for fabricating mouthguards, in a collage of funky colours.

Ref - FG



GLITTER GUARD BLANKS

A soft 4.00mm EVA blank ideal for fabricating mouthguards, in a collage of glittery colours.

Ref - GG



CLEAR RIGID: RETAINER SPLINT MATERIAL

This material is ideal for Retainers, Splints, Surgical Trays, Bite Raisers etc. It has excellent forming qualities, acrylic will bond to it and has a life span in excess of 6 months under normal circumstances.

Available in :	127mm Square	120mm Round	125mm Round
0.50mm Pk of 10 Blanks	VT5	VT5-120	VT5-125
0.50mm Pk of 100 Blanks	VT5-L	VT5-L-120	VT5-L-125
0.75mm Pk of 10 Blanks	VT75	VT75-120	VT75-125
0.75mm Pk of 100 Blanks	VT75-L	VT75-L-120	VT75-L-125
1.00mm Pk of 10 Blanks	VT1	VT1-120	VT1-125
1.00mm Pk of 100 Blanks	VT1-L	VT1-L-120	VT1-L-125
1.50mm Pk of 10 Blanks	VT15	VT15-120	VT15-125
1.50mm Pk of 100 Blanks	VT15-L	VT15-L-120	VT15-L-125
2.00mm Pk of 5 Blanks	VC2	VC2-120	VC2-125
3.00mm Pk of 5 Blanks	VC3	-	-



Note this material is hydroscopic and therefore can absorb moisture if not stored correctly, we recommend you buy no more than three months stock at any one time.

PROFORM DUAL LAMINATE

The dual laminate is a unique combination of two materials formed into a single laminated sheet. It was made especially for use in bruxing cases. The soft liner assures patient comfort, the hard surface assures long-lasting wear. Also, acrylic will bond to the hard surface to allow the dental technician to modify the occlusal surface.

Available in :	127mm Square	120mm Round	125mm Round
3.00mm sold per blank	DLHS3	DLHS3-120	DLHS3-125
5.00mm sold per blank	DLHS5	DLHS5-120	DLHS5-125



AQUEOS SPORT MOUTHGUARD DISINFECTANT

Safe to use on all mouthwear and containers without causing distortion or degradation. Kills 99.999% bacteria, fungi & viruses in less than 30 seconds.

60ml Ref: 48SD



PROTHO-CLEAN

Synthetic cleaner for high lustre finish, ideal for dentures, orthodontic and mouthguard appliances etc. Just spray on, soak for a second and wipe off with a soft towel. Removes polish and monomers leaving them with a high sheen, improves presentation.

75ml spray Ref: 04PCS





SILICONE
RELEASE SPRAY

Silicone spray eliminates sticking problems, we recommend its use on plaster models, units clamping frame and vacuum bed, it also prevents bonding between certain types of plastic.

Ref: 50CSS



FREEZE SPRAY

Use this spray immediately after vacuum forming and whilst the vacuum is still switched on, so the spray cools and contracts the plastic rapidly to ensure you achieve an accurate fit to the model. The spray is non flammable.

Ref: 49CFS



KEYSTONE
MACHINE III

Dual chamber vacuum forming unit, economically priced, suitable for most vacuum forming applications.

Ref: 05VAU



TRI-CUTTER

3 sided pyramid shape

Each

Ref: 40TC



TWIST
DRILL BUR

Drill bit shape

Ref: 07JLC



LARGE
DIARIT

For retainers

Ref: 33DTL



ARROW GEL

Burns out without residue

Green 3g Ref: 59AGG

Red 3g Ref: 59AGR

Centri-fuse Kit	Ref: CFSK
Refill pack of 25 centri-fuse strips	Ref: CFRP
Primer 10 mls	Ref: CFP
Flameless Torch	Ref: MHG
Pack of 12 Goat Hair Mounted Brushes	Ref: GH21
Scissors	Ref: FC
Decal Applicator Brush	Ref: CFHB
Pack of 5 A4 Inkjet Water Transfer Paper	Ref: CFIWTP
Pack of 5 A4 Sealing Paper for Inkjet Paper	Ref: CFSP
Pack of 5 A4 Laser Water Transfer Paper	Ref: CFLPP
Pack of 12 Pre-Printed Decals quote decal ref	Ref: CFPPD



SCOTCH BRITE DISCS

These discs are ideal for polishing and finishing mouthguards and other thermoplastic formings.

Available in :	Pk of 10	Pk of 100
Brown Coarse	HSBC	HSBCL
Grey Medium	HSBM	HSBML
Purple Fine	HSBF	HSBFL
Orange Extra Fine	HSBXF	HSXBFL
White Non Abrasive	HSBS	HSBSL



GREY SILICONE POLISHERS

Ideal for finishing and polishing most thermoplastics, supplied mounted.

Available in :	Bullet Shape	Large Domed
Pack of 5 Mounted Polishers	03BUP8	03BUP9
Pack of 100 Mounted Polishers	03BUP8L	03BUP9L



BLOCK-OUT RESIN

A low viscosity, blue pigmented light curing block-out resin for blocking out undercuts or creating reservoirs for the bleaching technique.

Pack of 4 Syringes	Ref: 01LCB0
Pack of 20 Syringes	Ref: 01LCB020



CENTRI-FUSE

An innovative product that makes the application of decals to mouthguards quicker, easier and more cost effective. The centri-fuse system eliminates the necessity of a second forming. The material is user friendly and adjustments are made in a matter of a few minutes with conventional tools.

CENTRI AFTERCARE

JEFE SILICONE DENTURE BRUSH

Unlike a regular brush, the soft silicone bristles are specifically made for dentures and flexible dentures. The Jefe is non abrasive and will help keep your denture in excellent condition.

Each Ref: 01STB

OLIVAFIX

A denture fixative cream like no other. What is so special about this product? It's simple. It is the first denture fixative cream on the market that contains NO Zinc & NO Petrolatum (more commonly referred to as "Vaseline"). The Zinc and Petrolatum have been replaced by almost 35% Olive Oil. Combining this with a strong, long lasting hold for up to 24 hrs and the fact that Olivafix contains no added

NITRADINE SHINE

With its "Anti-biofilm" formula, NitrAdine SHINE has been especially formulated to clean, whiten, refresh

NITRADINE

Most removable dental appliances (dentures) are made from porous materials in which opportunistic and pathogenic micro-organisms (fungi, bacteria and viruses) infiltrate, attach and form bio-film, a 3-dimensional protective matrix in which they build-up. Especially, Candida albicans penetrates deeply the pores and cracks of the dental appliance. Also other pathological micro-organisms such as Pseudomonas St mutans, St aureus (inc MRSA), aeruginosa and Herpes simplex can attach and colonise the appliance which leads to painful irritation of the oral cavity (denture stomatitis).

PROFESSIONAL DENTURE
CLEANING TABLETS

Our professional effervescent denture cleanser tablets are simple, efficient and easy to use. The combination of the active ingredients removes plaque, impurities, food remains, nicotine, tea, coffee and other stains.

Per 32 tablets Ref: 01PNT

TCS FLEXIBLE DENTURE CLEANER

The perfect companion to your new flexible denture, the TCS flexible cleaner - CLEANS, DISINFECTS & REFRESHES!

10 months' supply Ref: 01TCSDC

CENTRI 5

Complete after care solution for any acrylic denture. The kit has everything the patient needs for 5 weeks.

- Centri toilet bag
- Denture box

- Silicone toothbrush
- NitrAshine
- Olivafix
- Cleaning tabs (box)
- NitrAdine tabs (10)
- Booklet

CENTRI 8

Complete after care solution for any flexible denture. The kit has everything the patient needs for 8 weeks.

- Centri toilet bag
- Denture box

- Silicone toothbrush
- NitrAshine
- Olivafix
- TCS Cleaner (2)
- NitrAdine tabs (box)
- Booklet



PACKAGING

Free silver print on all boxes (subject to minimum order qty's)

HIGH GLOSS ORTHO BOXES PK 10



DENTURE BOXES PK 10



CROWN BOXES

1" crown box
2 Foams - Per 100

Ref: 01CBS

2" crown box
2 Foams - per 100

Ref: 01CBM



CROWN BOTTLES

16x57mm
Per 1,000

Ref: 01CDBS

26x61mm
Per 500

Ref: 01CDBL



INKO PACK

Pk 25

Ref: 39MST



APPLIANCE BOXES

Small

Pink Blue
Per 10 Per 10

Ref: 01VAC Ref: 01VACB

Large

Pink Blue White
Per 10 Per 10 Per 10

Ref: 01ABP Ref: 01ABB Ref: 01ABW



SLIMLINE BOX WHITE

Per 10

Ref: 01SLBW



DR DENTI BOXES

With brush and mirror. Each

Ref: 07DDDC



DENT SPA

Vibrating denture box. Each

Ref: 01DSC



TOOTH FAIRY BOX

Glow in the dark. Per 10

Ref: 39GITDB

CENTRI

COLACRYL



System One: Cold cure, chair side, soft reline tissue conditioner - 150g / 100ml kit.

Ref:01CSRS



System Two: Cold cure, chair side, hard reline for permanent relines - 150g / 100ml kit.

Ref:01CHRS

OPACRYL ORTHODONTIC ACRYLIC

Specifically designed and produced for the spray technique, this effective self-curing acrylic does not run or slump and represents great value for money. Available in Rose or Clear and 9 feature colours in 125ml: violet, orange, turquoise, black, blue, rubine, yellow, green and red.



1kg Polymer	Ref: 010P1
3kg Polymer	Ref: 010P3
500ml Monomer	Ref: 010L50
1ltr Monomer	Ref: 010L1
2.5ltr Monomer	Ref: 010L2
Std Pk (1kg/500ml)	Ref: 010SP
125ml Coloured Liquid	Ref: 010CL

CENTRIBASE ACRYLIC DENTINE'S



Available in 12 Acrylite shades and some Vita™ shades these heat cure and cold cure dentines are a quick and easy solution to producing temporary teeth.

100g HC or CC

CENTRIBASE SPECIAL TRAY



1kg Polymer	Ref: 01STP1
3kg Polymer	Ref: 01STP3
500ml Monomer	Ref: 01STLS
1ltr Monomer	Ref: 01STL2
Std Pk (1kg/500ml)	Ref: 01STMS
Bulk Pk (3kg/2ltr)	Ref: 01STML

CENTRIBASE DENTURE BASE



Available as standard (heat cure) or centribase X with an "extended dough time" - the material is available in 4 shades. Clear (glass-std only), light pink veined, translucent & natural veined. Powder can be supplied un-veined upon request.

454g Polymer	Ref: 01FCP1
3kg Polymer	Ref: 01FCP3
4 x 3kg Polymer	Ref: 01FCP12
8 x 3kg Polymer	Ref: 01FCP24
50kg Polymer	Ref: 01FCP50
1ltr Monomer	Ref: 01FCL1
2.5ltr Monomer	Ref: 01FCL2

Kits:	
Laboratory Kit (3kg / 1ltr)	Ref: 01FCPL
Intro Kit (3 x 3kg & 2 x 2.5ltr)	Ref: 01FCBP
Bulk Kit (8 x 3kg & 4 x 2.5ltr)	Ref: 01FCBP

Cold cure is also available

CENTRIBASE POURABLE



A self-curing pourable acrylic suitable for fabricating precision dentures, relines and repairs. Mixes quickly to form a perfect consistency for pouring, with a good working time. Finishes and polishes to a high gloss with a high colour stability. Available in three shades light pink, translucent and natural in either veined or unveined.

454g Powder	Ref:01PCC45
500ml Liquid	Ref:01PCCL

DIAMOND D HI IMPACT



Diamond D Ultra Hi Impact acrylic has exceptional flexural strength that comes in three shades. Original, Light and Chroma shades are available in 3 pack sizes offering great value for money. Printed denture boxes and bags offer the complete solution.

Original 450g/236ml	Ref: 39DDHOSK
Original 2.25kg / 946ml	Ref: 39DDHOMK
Original 11.25kg p/only	Ref: 39DDHOLP
Light 450g/236ml	Ref: 39DDHLSK
Light 2.25kg / 946ml	Ref: 39DDHLMK
Light 11.25kg p/only	Ref: 39DDHLPL
Chroma 450g/236ml	Ref: 39DDHCSK
Chroma 2.25kg / 946ml	Ref: 39DDHCMK
Chroma 11.25kg p/only	Ref: 39DDHCLP

236ml Monomer	Ref: 39DDHL2
946ml Monomer	Ref: 39DDHL9

Cold cure is also available



TCS SUCTION CUP

TCS Suction Cup is a soft silicone rubber that can be used for maxillary and mandibular restorations. Existing dentures can easily be converted into TCS Suction Cup dentures. TCS Suction Cup utilises a series of tiny suction cups that adhere to the oral tissue without irritation. The silicone material is also an excellent permanent soft relene.

TCS Suction Cup Primer 5ml	Ref: 01TCSpsc
TCS Suction Cup Bur	Ref: 01TCSscb
TCS Suction Cup Mandrel	Ref: 01TCSscm
TCS Suction Cup Silicone Paste	Ref: 01TCSscp



TECNOFLASK

The Tecnoflask is a microwaveable flask that helps create acrylic dentures quicker and more efficiently. Suitable for all Hi Impact and standard acrylic denture base. The Tecnoflask will give your laboratory a cleaner environment and make the process more enjoyable.

The entire processing time can be completed in 10 minutes which represents in some cases a time saving of over 2 ½ hours. The starter kit also comes with a press board to prevent damage to the flask when using it with a clamp.

Each Ref:39TFIK



ITSOCLEAR

High Quality Crystal Clear Clasp, ItsoClear provides you with the ability to produce and repair a clasp in under 4 minutes with a cosmetic appeal unlike any other clasp.

- Saves bench time producing a higher quality dental aesthetic clasp that your customers will appreciate while increasing your profit.
- Crystal clear, will not lose its clarity
 - Non-abrasive
 - Easily adjustable
- No injection machines or hand-held injectors required
 - Perfect for new partial dentures as well as repairs and replacing metal clasps

Per 6 Universal	Ref:01EZCU	Per 6 Molar	Ref:01EZCM
-----------------	------------	-------------	------------



CENTRI COMPRESS

With unique spring system that significantly increases the longevity of the unit

Each 2 Flask	Ref: 01FC2F
Each 3 Flask	Ref: 01FC3F



CENTRI FLASKS

Each Flask	Ref: 05MSLF
Each Flask	Ref: 23F89BS



CENTRI HYDRAULIC PRESS

Each 3 Flask	Ref: 55CHP
--------------	------------

ALPHABETICAL PRODUCT LISTING

A Acid Etch..... 10 Aegis Crown and Bridge Resin Cement 32 Aegis Liner 33 Aegis Pit and Fissure Sealant..... 29 Aegis V 34 Apavia AVT Max 8 Apavia EcoMax 8 Apavia HydroMax..... 9 Apavia PurMax 9 Apavia UltraMax..... 9 Appliance Boxes 42 Aqueous Sport Mouth Disinfectant 39 Arrow Gel..... 40	G Gelato APF Foam Fluoride..... 26 Gelato APF Fluoride Gel 26 Gelato Home-care Stannous Fluoride Gel & Vivid Neutral Sodium Fluoride Gel..... 27 Gelato Topical Anaesthetic Gel 11 Gelato Benzocaine Anaesthetic Metered Spray 11 Gelato Perio Maintenance Rinse..... 26 Gelato Proply Paste 28 Glaze 15, 25, 31 Glitter Guard Blanks..... 39 Gorilla Seal 25, 31, 34 Grey silicone Polishers..... 40 G.T.C. 19	Professional Denture Cleaning Tablets..... 41 Proform Dual Laminate 39 Proform Nightguard..... 36 ProplyBrite..... 28 Proply Powder..... 28 Prothro-clean..... 39
B Big Bite Tray..... 18 Bleach Tray Material..... 36 Block-out Resin..... 40 Bo Boxes..... 25 Bonding Liquid (under New Truliner & Trusoft)13, 14 Bookkeeping Forms 7 BoTrays..... 20	H Handles 7 High Gloss Ortho Boxes pk 10..... 42 Hinged EZ Fluoride Trays..... 27 Hi-Spot..... 12 Hydrox 33	R Re'CORD 18 Rimseal 22
C Centri 5..... 41 Centri 8..... 41 Centribase Acrylic Dentine's..... 43 Centribase Denture Base..... 43 Centribase Pourable 43 Centribase Special Tray..... 43 Centri Compress 44 Centri Flasks..... 44 Centri Hydraulic Press 44 Clear HD for Nightguards and Fluoride Trays 36 Clear Soft Mouthguard Blanks 38 Clear Rigid: Retainer - Splint Material 39 Colacryl 43 ComfortCaine 10 Cooley Pic..... 5 Copaliner 33 Coping Material..... 36 Cora-Caine..... 10, 12 Crown Bottles..... 42 Crown Boxes..... 42	I Impact 1500..... 11 Impact 2000..... 12 Inko Pack..... 42 ItSoClear 44	S Sapphire 22 Scotch Brite Discs 40 Silicone Release (under DentuSil) 13 Silicone Release Spray 40 Simplastic Dispenser 16 Simplastic Hazard Bag 16 Simplastic Patient Drape 16 Simplastic Sleeve..... 16 Simplastic Tray Cover 16 Simplastic Uni-Cover 16 Softone..... 15 Slimline Box White..... 42 Small Diarit..... 40 Spatulas..... 7 Special Designs Mouthguard Blanks..... 38 Stainless Steel Trays 21 Superbite..... 18 SuperEBA 32 Supergel..... 22 Supergel Fresh 22 Supergel Magic 22 Superpaste..... 23 Supersil..... 23 Syringes 7
D Dental Necessities..... 7 Dent Spa 42 DentuSil 13 Denture Boxes 42 Diamond D Hi Impact..... 43 Dr Denti Boxes..... 42 Duz-All 14	J Jefe Silicone Denture Brush..... 41	T TAC 17 TC Retractors 7 TCS Flexible Denture Cleaning..... 41 TCS Suction Cup..... 44 Tecnoflask 44 Tongue Scrapers..... 29 Tooth Fairy Box 42 Tray Aways..... 20 Tray and Base Plate Material..... 36 Trim 30 TrimII 30 Trim Plus 30 Tri-cutter 40 TruRepair 15 TruSoft 14 Twist Drill Bur..... 40
E Edentulous Tray Aways 20 Esteem 34 Evacuators 5 E-Z Squeeze Universal Dispenser..... 5	K Keystone Machine III..... 40	U Ultra Fast 15 Ultra Trim 31
F Fastray..... 19 Fastray LC 19 Flexo..... 6 Foam Backed Bleach Tray..... 36 Fortis Disposable Proply Angle 27 Freeze Spray 40 Funguard Blanks 39	L Light Liner..... 13 Large Diarit 40	V Vigilance 16 Visual Systems 7
G Gelato APF Foam Fluoride..... 26 Gelato APF Fluoride Gel 26 Gelato Home-care Stannous Fluoride Gel & Vivid Neutral Sodium Fluoride Gel..... 27 Gelato Topical Anaesthetic Gel 11 Gelato Benzocaine Anaesthetic Metered Spray 11 Gelato Perio Maintenance Rinse..... 26 Gelato Proply Paste 28 Glaze 15, 25, 31 Glitter Guard Blanks..... 39 Gorilla Seal 25, 31, 34 Grey silicone Polishers..... 40 G.T.C. 19	M Masque 6 Masterlink..... 32 Mixing Brushes 5 Mixing Cups 5 Mixing Pads..... 5 Mizzy Prevent Trays..... 27 "My Tooth" Boxes..... 6 Multi-Coloured Striped Mouthguard Blanks 38	W Wonder Wedges 7
H Handles 7 High Gloss Ortho Boxes pk 10..... 42 Hinged EZ Fluoride Trays..... 27 Hi-Spot..... 12 Hydrox 33	N Nite Bite 29 NitrAdine 41 NitrAdine Shine 41 Neocryl 24 Neocryl Splint 24 New Truliner 13 Novelty Keychains..... 6	Z Zinc Phosphate 33 Zeta..... 35
I Impact 1500..... 11 Impact 2000..... 12 Inko Pack..... 42 ItSoClear 44	O Olivafix..... 41 Opacryl Orthodontic Acrylic 43 Orange Solvent 19 Original Truliner 14 OrthoGlo 25	
P Pit & Fissure Sealant with Fluoride 29 Plastogum 23 Plastopaste 22 Pressure Pot 6		

CLINICAL & CDT CATALOGUE

CONTACT

WHW Plastics Ltd
Therm Road, Cleveland Street
Hull, HU8 7BF
United Kingdom

Tel: +44 (0) 1482 329154
Fax: +44 (0) 1482 217140

Email: sales@whwplastics.com
Web: www.whwplastics.com