

How many restorations can be milled per cutter/bur?

	Cutter					
	0.6 mm	1.0 mm	2.5 mm	0.3 mm	TI 1.0 mm	TI 2.0 mm
Ceramill Zolid	90			75		
Ceramill ZI	90			75		
Ceramill Wax	200			1200		
Ceramill PMMA/TEMP	45			/		
Ceramill TEMP Multilayer	15			/		
Ceramill Splintec	/	10	10	/		
Ceramill Sintron	200			240		
Ceramill PEEK	130			/		
Ceramill TI-Forms					25	25



Stated quantities can differ depending on indication, size of block and individual geometry!

How many restorations can be milled per cutter/bur?

	Bur (diamond, orange)			
	1.0 mm	1.4 mm	1.8 mm	0.4 mm
VITA Enamic	100	100	30	10
VITA Mark II	50	50	30	10
VITA TriLuxe forte	50	50	30	10
e.max CAD	15	15	12	10
VITA Suprinity	25	25	25	10
3M Lava Ultimate	40	40	40	10



Stated quantities can differ depending on indication, size of block and individual geometry!

How many restorations can be milled per cutter/bur?

	Roto DMB DC			
	0,6 mm	1,0 mm	2,5 mm	0,3 mm
Yamakin KZR-CAD HR 2	60	80	80	
GC Cerasmart	60	80	80	
Shofu Block HC	80	100	100	80



Stated quantities can differ depending on indication, size of block and individual geometry!