



The Last Christmas Ticket Competition is bowing out in spectacular style as we will be giving away a holiday for two, to Bali, Indonesia for 7 nights. Staying at a 4 star hotel with flights included this really is an opportunity of a lifetime.

The most popular of the Indonesian Archipelago's and adorned with white sandy beaches, Bali is most known for its forested volcanic mountains, iconic rice paddies and coral reefs. Leave the wax and plaster behind and just relax.

entries you can receive.

The Last Christmas Ticket Competition draw will take place on Wednesday December 20th @ 5pm. All entries up until this time will appear on your invoice. For full terms and conditions please look at



Offers will be valid from Dec 1st 2017 until Jan 31st 2018 (The Last Christmas Ticket Competition entries will cease as of purchases from December 21st onwards)





## **SILENT** powerCAM EC

Bag-less, high-performance extraction unit for CAM systems.

- O Low-maintenance brushless motor with a guaranteed 5,000 hours of operation
- O Suitable for CAD-CAM systems requiring high suction performance due to a powerful motor with a maximum flow rate of 4000l/min
- Sophisticated bidirectional communication with the CAM system via PLC control
- O No dust bags to change due to bag-free fine filter technology, including automatic filter cleaning
- O Very low operating F noise of 54.3 dB(A)

## **SILENT EC2**

Bag-less, highperformance automatic extraction unit for two benches

R2939000003

- O Low-maintenance brushless motor with a guaranteed 5,000 hours of operation
- Consistently high suction thanks to the intelligent power control
- No dust bags to change due to bagfree fine filter technology, including automatic filter cleaning
- Very low operating noise of 55 dB(A)
- Tuneable autoextraction



## **High-performance** technology for professionals

#### **Optimised performance**

Suitable for the most difficult grinding, separating and milling operations, the Q Profi manages them all with ease. With an optimisable speed range of between 1,000 - 50,000 rpm (knee version) and 8.0 Ncm, the Q Profi represents power, efficiency and style.

The soft grip, precisely balanced handpiece creates the benchmark in comfort and control The shape of the handpiece and the integrated sealing system provides excellent protection against pollution and premature wear

99SCH 9340/5 Q Profi Knee Control Complete 99SCH 9350/5 Q Profi Foot Control Complete





Efficient compact extraction unit for CAM systems.

- Sophisticated bidirectional communication with the CAM system via SPS interface
- O Especially developed for the CAM market
- No dust bags to change due to bag-free fine filter technology, including mechanical filter cleaning
- O Low-maintenance motor with a guaranteed 1,000 hours of operation







Compact extraction for either one appliance or a single workplace.



- Two suction levels
- O No dust bags to change due to bagfree fine filter technology, including mechanical filter cleaning
- O Low-maintenance motor with a guaranteed 1,000 hours of operation
- O Tuneable autoextraction





8 · × / • •



Polyflex is a reversible hydrocolloid material to be used for duplicating any model in stone or phosphate investments. The material is an ideal hydrocolloid for fluid resin techniques and can also be used with ethyl-silicate investments.

With incredible stability it's the perfect partner for Lucitone® HIPA high impact pourable acrylic and your basic laboratory

needs









Used for the acrylic pour technique and/or to duplicate models. Manufactured in light-

## Lascod Xilgum 70

Used for the fast reproductions of gum morphology

- O Scannable
- Short setting time
- O High detail definition (20µm)
- O High dimensional stability
- O Suitable for direct and in-direct technique
- O No shrinkage or deformation
- O Easy to finish
- Natural colour
- O Working time @23°C -2 minutes
- O Setting time @ 23° 10 minutes
- O Dimensional shrinkage @ 24 hours 0.02%

## Lascod Ghenesyl 92

- O Intra-oral impression taking A-Silicone putty
- Next generation polyvinylsiloxane giving excellent definition
- Incredibly thixotropic and hydrophilic
- O Fast setting time in mouth
- Tear resistant
- O Elastic
- Biocompatible
- O Adaptable to the anatomical shapes reproducing perfectly integral margins

## **Lascod TIME OUT**

- Liquid retardant able to modify the working time of plasters/stones
- O TIME OUT does not affect the plaster characteristics
- O Extremely useful for casting multiple models at the same time
- O Every drop of TIME OUT added to the water when using synthetic gypsums increases the working time by approximately 20 seconds...



66LTO

Xilgum

Goosperia Goodperie di Goodperie Propus Propus Liberatione Per Midger Per Midger Midger User Midger Midger User Midger User Midger Midger User Midger User Midger

SYSTEM



lilgun



## **Verty High Security Fluid 75**

- O Silicone insulating material
- O Protects the acrylic or porcelain teeth during investing
- O Easily applied in one single coat onto the teeth and wax
- O Forms a thin, very hard insulating layer
- O Mixing time 45 sec
- O Working time 4 min
- O Setting time 15 min





Snoring is caused from a blockage of the tension in the muscles is released which causes the throat to narrow. The respiratory air is forced through the narrowed throat channel and the soft palate begins to vibrate, producing the typical snoring noises. The tongue can then fall back into the throat and this might lead to a complete blockage which can result in the periodical interruption of breathing (apnoea).

DocSnoreNix® is a newly developed and patented SNX-protrusion system according to Dr. Schmitt-Bylandt, which brings the patient's lower jaw into an

**DocSnoreNix® Starter Set Contents** 

O 1 x patient set of protrusion springs

O Thermoforming blanks (2 x

Kombiplast 3 mm)

Protrusion calibre

transparent

ml bottle

Accessories

Manual

O Lightdon Gel, 20 ml

O Lightdon Bonding, 10

O Ultratrimm medium

Pair of pliers

individually determined protruded position during respiratory system during sleep. While sleeping, the sleep. By doing so the throat is widened and the typical danger of blocking is decisively reduced.

> Unique, super-elastic connecting elements between the upper and lower jaw allow any kind of movement of the jaw. The strain on the mandibular joint is significantly reduced using the delicately

designed "spring elements". This creates a comfortable non-intrusive appliance which means a greater acceptance and use by the patient.



- O 1 x patient set of protrusion springs
- O Thermoforming blanks (2 x Kombiplast 3 mm)
- O Lightdon Gel, 20 ml transparent
- O Lightdon Bonding,



Book now for your one day **DocSnoreNix®** course on either March 12th, 14th or 16th 2018. Become a certified DocSnoreNix manufacturer for only £325 +VAT with 8 hours verifiable CPD





We are pleased to announce the start of our trade in program on the purchase of certain new equipment lines. When you see this logo. We will consider all exchanges on the same type of product.

#### Monotrac V2

When it comes to model, die and articulation, "the cost, level of efficiency and accuracy" are the constant questions. Monotrac is the answer. Offering more options, accuracy and efficiency than any other model system out there.

The integrated articulator, with pinless cone matrix model base, allows a single pour "cast-in-place" set-up that is simple and fast while eliminating labour intensive pinning and basing procedures.

Monotrac dies are easily removed and reseated with solid stable precision. Even lower anteriors, the thinnest of dies, resist lateral sway. Other key features of the Monotrac are the Flex Arm Hinge allowing natural excursive equilibration. The Plug-in adjustable vertical stops are easy to place and eliminate the need for a metal articulator set up. One system covers it all from refractory models to sectional implant full arches, it can do it all.



- O 18 V2 Straight Quad Bases
- O 16 V2 Anterior Bases
- O 5 V2 Full Arch Bases
- O 12 V2 Radius Quad Bases
- O 28 V2 Hinge Assembly (sets)
- 6 Plug in adjustable stops
- O 1 Hammer
- O Complete step by step instructions for all 4 setup options

#### **Monotrac Sampler Kit**

- O 16 V2 Straight Quad Bases
- O 4 V2 Anterior Bases
- 0 4 V2 Full Arch Bases
- O 13 V2 Hinge Assembly (sets)
- 12 Plug in adjustable stops
- O Complete step by step instructions



53MTV2SK









#### Renfert EASY view 3D

"Making work easy" for dental technicians is our mission at Renfert. With our new EASY view 3D Dental Viewer, we have developed an innovative video microscope with a 3D monitor that achieves this objective from a whole range of different perspectives. EASY view 3D helps users maintain an ergonomic posture while enabling fast, direct communication between technicians, the lab and the dentist's practice as well. Renfert's EASY view 3D Dental Viewer is not only ideally suited for viewing objects, for routine processing and quality control, but equally useful for educational, training and documentation purposes.

**Precise and reliable:** 15x and 20x magnification in 2D or 3D mode

**Ultra-high definition:** Image reproduction on a 21.5-inch full HD 3D monitor

**All details at a glance:** The passive polarization glasses enable precise 3D viewing and are well suited for persons wearing glasses

The EASY view 3D's three-dimensional mode enables easy and natural hand-eye coordination. Its three-dimensional spatial representation is perfect for critically examining work objects. This helps streamline your entire production process and makes many routine tasks, like finishing and performing quality checks,

much more precise. Depending on the indication, a simple mouse click lets you toggle between the 3D and 2D modes. The 2D mode is ideal for checking and viewing objects.

Material-optimized image reproduction through 4 different object modes. Toggle through the different modes by mouse click.

# C MODE (CERAMIC)→ M MODE (METAL)→ S MODE (METAL)→ N MODE (NEUTRAL)→

A dental technician's work assumes a new dimension with the EASY view 3D Dental Viewer. The ultimate workplace ergonomics enable a sustained natural and healthy posture. This enhances the users' well-being and performance. The flexible boom arm can be adjusted to any work situation and object. The monitor has been especially developed for system compatibility and ensures sharp high-definition image reproduction. Positioned behind the object, the monitor is ergonomically and ideally located for an unobstructed view of the work object.

The system is conveniently operated by mouse and footswitch, thereby ensuring a relaxed and healthy working posture.

The easy-to-use and intuitive functionality of the EASY view 3D Dental Viewer helps you document work, share information and data with your team as well as the dentist's practice. Catch the monitor image by footswitch and store the screenshot on the USB stick or send it via LAN connection to your PC. Images and video sequences are stored in 2D mode.

CALL NOW FOR A
DEMONSTRATION OF
THIS REVOLUTIONARY
CONCEPT



# TOERING FIBERS









**CENTRI BLADE** 







We will be closing on Thursday December 21st 2017 until Tuesday January 2nd 2018, orders that need to be shipped out on the 21st need to be received by us before 3pm.

From everyone at WHW we wish you a

Merry Christmas, Happy Holiday and a joyous New Year!

#### **WHW Plastics Ltd**

Therm Road, Hull, East Yorkshire, HU8 7BF United Kingdom tel +44 (0)1482 329 154 • fax +44 (0)1482 217140

• email: sales@whwplastics.com • website: www.whwplastics.com

Offers valid from December 1st 2017 until January 31st 2018















